TABLE DES MATIÈRES

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TABLE OF CONTENTS

	INTRODUCTION	140
I	WARNINGS	140-143
	CONTRAINDICATIONS	140-141
	SAFETY MEASURES	141-143
II	PRESENTATION.....	144-147
	RECEIPT OF EQUIPMENT AND ACCESSORIES.....	144
	WARRANTY.....	145
	MAINTENANCE	145
	STORAGE AND TRANSPORTATION CONDITIONS	145
	USE CONDITIONS.....	145
	DISPOSAL	146
	STANDARDS.....	146
	PATENTS	146
	STANDARDIZED SYMBOLS	146
	TECHNICAL CHARACTERISTICS.....	147
III	HOW DOES ELECTROSTIMULATION WORK?	148-149
IV	USAGE GUIDELINES	150-151
	PLACEMENT OF ELECTRODES	150
	STIMULATION POSITIONS.....	151
	ADJUSTING STIMULATION ENERGIES.....	151
V	HOW IT WORKS.....	152-160
	DESCRIPTION OF THE DEVICE	152
	BATTERY INSERTION.....	153
	CONNECTIONS	153
	PRELIMINARY SETTINGS	153-154
	SELECTING A CATEGORY.....	154
	ADJUSTING STIMULATION ENERGIES.....	157
	PROGRAM PROGRESSION.....	157-158
	END OF PROGRAM.....	158
	BATTERY LEVEL AND RECHARGING	158
	PROBLEMS AND SOLUTIONS.....	159-160
VI	PERSONAL TRAINER GOALS.....	161-163
VII	WELLNESS PROGRAMS	164-168
VIII	EMC TABLES.....	169-172

INTRODUCTION

Please read this manually carefully prior to using your Compex Fit stimulator. You are strongly advised to read Chapter 1 of this manual.

The Compex Fit is a stimulator intended for muscular training and pain relief. Any person excepted those mentioned in the Chapter 1 "Warnings" can use the Compex Fit.

EN

I WARNINGS

CONTRAINDICATIONS

Important contraindications

- Cardiac stimulator (pacemaker)
- Epilepsy
- Pregnancy (do not use on abdominal region)
- Serious arterial circulation disorders in lower limbs
- Abdominal or inguinal hernia

Precautions when using the Compex

- After trauma or a recent operation (less than 6 months)
- Muscular atrophy
- Persistent pain
- Need for muscular rehabilitation

Osteosynthesis equipment

The presence of osteosynthesis equipment (metallic equipment in contact with the bone: pins, screws, plates, prostheses, etc.) is not a contraindication for the use of Compex programs. The electrical currents of the Compex are specially designed to have no harmful effect on osteosynthesis equipment.

CONTRAINDICATIONS

Important!

- Do not use Compex stimulator programs if you have sensitivity problems.
- Never use the Compex for prolonged periods without medical advice.
- Consult your doctor if you are in any doubt whatsoever.
- Read this manual carefully, in particular Chapters 6 and 7, which provide information concerning the effects and indications of each stimulation program

EN

SAFETY MEASURES

What you should not do with the Compex

- Do not use the Compex in water or in a humid atmosphere (sauna, hydrotherapy, etc.).
- Do not use the Compex in an oxygen-rich atmosphere.
- Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, very highly strung people may have a vagal reaction. This is of psychological origin and is connected with a fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. This vagal reaction produces a feeling of weakness with a fainting tendency, slowing down of the heart rate and a reduction in arterial blood pressure. If this does occur, all that is required is to stop the stimulation and lie down with the legs raised for the time it takes for the feeling of weakness to disappear (5 to 10 minutes).
- Never allow movement resulting from muscular contraction during a stimulation session. You should always stimulate isometrically; this means that the extremities of the limb in which a muscle is being stimulated must be firmly fixed, so as to prevent the movement that results from contraction.
- Do not use the Compex if you are connected to a high-frequency surgical instrument as this could cause skin irritation or burns under the electrodes.
- Do not use the Compex within X meter (see EMC table) of short wave or microwave devices as this could alter the currents generated by the stimulator. If you are in any doubt as to the use of the Compex close to another medical device, seek advice from the manufacturer of the latter or from your doctor.

- Do not use the Compex in areas where unprotected devices are used to emit electromagnetic radiation. Portable communications equipment can interfere with electrical medical equipment.
- Only use stimulation cables supplied by Compex.
- Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first.
- Never connect stimulation cables to an external power supply. There is a risk of electric shock.
- Never use a battery other than the one supplied by Compex.
- Never recharge the stimulator without first disconnecting the stimulation cables.
- Always use the charger provided by Compex to recharge the batteries.
- Never use the Compex or the charger if it is damaged (case, cables, etc.) or if the battery compartment is open. There is a risk of electric discharge.
- Disconnect the charger immediately if the Compex "bleeps" continuously, if there is abnormal heating or smell, or if smoke comes from the charger or the Compex.
- Do not recharge the battery in a confined space (carrying case, etc.). There is a risk of fire or electric shock.
- Keep the Compex and its accessories out of reach of children.
- Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the Compex, the battery compartment and the charger.
- Sudden temperature changes can cause condensation to build up inside the stimulator. Use the device only once it has reached ambient temperature.
- Do not use the Compex while driving or operating machinery.
- Do not use the stimulator at altitudes of over 3,000 metres.

Where never to apply the electrodes

- In the vicinity of the head.
- Contralaterally, i.e.: do not use two poles connected to the same channel on opposite sides of the body.
- On or near skin lesions of any kind (wounds, swelling, burns, irritation, eczema, etc.).


EN

Precautions when using the electrodes

- Only use electrodes supplied by Compex. Other electrodes may have electrical properties that are unsuitable for the Compex stimulator.
- Always turn off the stimulator before moving or removing any electrodes during a session.
- Do not place the electrodes in water.
- Do not apply solvents of any kind to the electrodes.
- For best results, wash and clean the skin of any oil and dry it before attaching the electrodes.
- Attach the electrodes in such a way that their entire surface is in contact with the skin.
- For obvious reasons of hygiene, each user must have his/her own electrode set. Do not use the same electrodes on different people.
- Never use a set of adhesive electrodes for more than 15 sessions as their bonding power deteriorates over time and optimal contact is very important for both user comfort and efficacy.
- Some people with very sensitive skin may experience redness under the electrodes after a session. Generally, this redness is completely harmless and disappears after 10 to 20 minutes. Never start another stimulation session in the same area, however, if the redness is still visible.

II PRESENTATION

EN

Receipt of equipment and accessories

Your kit contains :

- A** 1 stimulator (1525000)
- B** 1 charger (68302x)
- C** 1 set of white Snap connectors for electrodes, colour-coded (blue, green, yellow, red) (601136)
- D** 2 bags of small electrodes (5x5 cm) (6260760)
- E** 2 bags of large electrodes (5x10 cm) (6260770)
- F** 2 user's manuals (45251xx)
- G** 1 belt clip (949000)
- H** 1 travel pouch (6680036)
- I** 1 quick start guide (4525060)
- J** 1 battery (941210)

Warranty

See the enclosed leaflet.

Maintenance

Do not sterilise the stimulator. Clean using a soft cloth and an alcohol-based, solvent-free cleaning product. Use only a minimum amount of liquid when cleaning the Compex. Do not attempt to repair the stimulator or any of its accessories. Never dismantle the Compex or the charger containing high voltage parts because of risk of electric shock.

Chattanooga Group declines all responsibility for any damages or consequences resulting from unauthorized attempts to open, modify, or repair the stimulator. This may only be done by persons or repair services authorized by Chattanooga Group. Your stimulator does not require calibration.

All units are systematically tested and validated prior to distribution. Their characteristics do not vary under normal conditions. Nonetheless, as the Compex is a high-quality electrical instrument, its lifespan depends on the use that is made of it and the care and maintenance it receives during its lifetime. If your stimulator contains parts that seem worn or defective, please contact your nearest Chattanooga Group customer service centre regarding an upgrade. Medical and health professionals are obliged to service the unit in compliance with relevant national laws and regulations. This involves testing performance and safety parameters at regular intervals.

Storage and transportation conditions

The Compex contains a battery. For this reason, it must be stored and transported in accordance with the following conditions:

Storage and transportation temperature: - 20°C to 45°C

Maximum relative humidity: 75 %

Atmospheric pressure: 700 hPa to 1060 hPa

Use conditions

Temperature: 0°C to 40°C

Relative humidity: 30 % to 75 %

Atmospheric pressure: 700 hPa to 1060 hPa

Do not use in areas where there is a risk of explosion.

Disposal

The main purpose of the 2002/96/EEC Directive is to prevent the creation of waste electrical and electronic equipment (WEEE) and to reduce the amount of waste to be disposed of by encouraging reuse, recycling, and other forms of recovery. The pictogram of a crossed-out wheelee bin means that the equipment may not be thrown out with the household refuse, but must be specially collected. The equipment has to be given to a suitable collection point for processing. By doing so, you will be contributing to the safeguarding of natural resources and health.

Batteries must be disposed of in compliance with relevant national regulatory requirements.

Standards

The Compex complies with current medical standards.

To guarantee your safety, the Compex has been designed, manufactured, and distributed in compliance with the requirements of European Directive 93/42/EC on medical devices. The Compex also complies with the standard on general safety regulations for electromedical equipment, IEC 60601- 1. It also complies with the IEC standard for electromagnetic compatibility, 60601-1-2 and the standard on particular safety requirements for nerve and muscle stimulators, IEC 60601-2-10.

Current international standards require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation).

The Compex also complies with Directive 2002/96/EEC on waste electrical and electronic equipment (WEEE).

Patents

The Compex incorporates several innovations with patents pending.

Standardized symbols

Caution: Read the user manual or operating instructions (symbol no. 0434 IEC 60878).

The Compex is a class II device with internal electric power and type BF applied parts (symbol no. 5333 IEC 60878).

Waste electrical and electronic equipment (WEEE) marking according to standard EN 50419.

The On/Off button is a multi-function button (symbol no. 5009 IEC 60878)

Functions

On/Off (two stable positions)
Waiting or preparation for part of the unit
Off (turns system off)

Technical characteristics

General points

941210 Battery

Nickel metal-hybrid (NiMH) rechargeable (4.8 V / \geq 1200 mA/h).

68302x Chargers

Only battery chargers bearing the following information can be used to recharge the batteries supplied with Compex stimulators:

Europe/ UK/ US

683020/ 683022/ 683026

Type TR503-02-E-133A03/ Type TR503-02-U-133A03/ Type TR503-02-A-133A03

Input 100-240 VAC / 47-63 Hz / 0.5 A max.

Output 9V / 400 mA / 6 W

601136 White Snap-connection stimulation cables

Device connector: 6-pin

Electrode connector: Female Snap

Length: 1,500 mm

Protection index

IPXO (IEC 60529)

Neurostimulation

All electrical specifications are given for an impedance of 500-1000 ohms per channel.
Channels: four independent and individually adjustable channels that are electrically isolated from each other and earthed.

Impulse shape: constant rectangular current with pulse compensation to eliminate any direct current component to prevent residual polarisation at skin level.

Maximum pulse intensity: 100 mA.

Pulse intensity increments: manual adjustment of stimulation intensity from 0 to 694 (energy) in minimum increments of 0.5 mA.

Pulse width: 60 to 400 μ s.

Maximum electrical charge per pulse: 80 microcoulombs ($2 \times 40 \mu$ C, compensated

Standard pulse ramp-up time: 3 μ s (20%-80% of maximum current).

Pulse frequency: 1 to 150 Hz.

III HOW DOES ELECTROSTIMULATION WORK?

The principle of electrostimulation is to stimulate nerve fibres by means of electrical impulses transmitted by electrodes. The electrical pulses generated by Compex stimulators are high quality pulses - offering safety, comfort and efficiency - which can stimulate different types of nerve fibres:

1. the motor nerves, to stimulate a muscular response. The quantity and the benefits obtained depend on the stimulation parameters and this is known as electro-muscular stimulation (EMS).
2. certain types of sensitive nerve fibres to obtain analgesic or pain-relieving effects.

1. Motor nerve stimulation (EMS)

In voluntary activity, the order for muscular work comes from the brain, which sends a command to the nerve fibres in the form of an electrical signal. This signal is then transmitted to the muscle fibres, which contract. The principle of electrostimulation accurately reproduces the process observed during a voluntary contraction. The stimulator sends an electrical current impulse to the nerve fibres, exciting them. This excitation is then transmitted to the muscle fibres causing a basic mechanical response (= muscular twitch). The latter constitutes the basic requirement for muscular contraction. This muscular response is completely identical to muscular work controlled by the brain. In other words, the muscle cannot distinguish whether the command comes from the brain or from the stimulator.

The parameters of the Compex programs (number of impulses per second, contraction time, rest time, total program time) subject the muscles to different types of work, depending on the muscle fibres. In fact, different types of muscle fibres may be distinguished according to their respective contraction speed: slow, intermediate and fast fibres. The fast fibres will obviously predominate in a sprinter, while a marathon runner will have more slow fibres.

With a good knowledge of human physiology and a perfect mastery of the stimulation parameters of the various programs, muscular work can be directed very precisely towards the desired goal (muscular reinforcement, increased blood flow, firming up, etc.).

2. Stimulation of the sensitive nerves

The electrical impulses can also excite the sensitive nerve fibres to obtain an analgesic or painrelieving effect. The stimulation of the tactile sensitive nerve fibres blocks the transmission of pain by the nervous system. The stimulation of another type of sensitive fibres creates an increase in the production of endorphins and, therefore, a reduction in pain. With pain relief programs, electrostimulation can be used to treat localised acute or chronic pain as well as muscle pain.

Attention : Do not use the pain relief programs for a long period without medical advice.

EN

Benefits of electrostimulation

Electrostimulation is a very effective way to make your muscles work:

- with significant improvement of different muscle qualities
- without cardio-vascular or mental fatigue
- with limited stress on the joints and tendons. Electrostimulation thus allows a greater quantity of work to be done by the muscles compared with voluntary activity.

To be effective, this work must involve the greatest possible number of muscle fibres. The number of fibres that are working depends on the stimulation energy. It will therefore be necessary to use the maximum tolerable energy. The user controls this aspect of stimulation. The higher the stimulation energy, the greater the number of muscle fibres that are working and, therefore, the more significant the progress achieved. To maximize your results, Compex recommends that you complement your electrostimulation sessions with other efforts, like:

- regular physical exercise
- proper and healthy nutrition
- balanced lifestyle

IV USAGE GUIDELINES

The usage guidelines presented in this section should be considered as general rules. For all programs, it is recommended that you read carefully the usage information and advice presented in chapters 6 and 7.

Placement of electrodes

For optimal results, use the electrode positions recommended by Compex. To do this, refer to the pictures and pictograms shown on the inside cover of the manual.

Each stimulation cable has two poles:

A positive pole (+) = red connection

A negative pole (-) = black connection

A different electrode must be connected to each pole.

NOTE: It is possible and normal to have an electrode arrangement that leaves one electrode connection free from a cable.

Depending on the characteristics of the current, efficacy can be optimised in certain programs by placing the electrode connected to the positive pole (red connection) "strategically".

For the objectives in the Personal Trainer category, which force muscles to contract, it is important to place the positive electrode on the motor point of the muscle.

It is crucial to choose the right size of electrodes (large or small) and correctly position these on the muscle group you want to stimulate so as to ensure the efficacy of the treatment. Therefore, always use the size of electrodes shown in the pictures. Unless you have other specific medical instructions, always follow the placement directions in the pictures. Where necessary, look for the best possible position by slowly moving the positive electrode over the muscle until you find the point that will produce the best contraction or the most comfort for you.

Compex disclaims all responsibility for consequences arising from electrodes placed in other positions.

EN

Stimulation positions

This position will vary depending on the position of the electrodes, the muscle group you wish to stimulate, and the program you are using. For the objectives in the Personal Trainer category, which involve powerful muscular contractions, the muscle should always be stimulated in an isometric fashion. You must therefore fix the extremities of your limbs securely. In this way, you provide maximum resistance to the movement and prevent any shortening of the muscle during the contraction, which could create cramp pains and serious stiffness after the session. For example, when stimulating the quadriceps, the user should be in a seated position with the ankles fixed with straps to prevent extension of the knees.

EN

For programs in the Wellness category, which do not induce powerful muscle contractions, but only muscle twitches, position yourself as comfortably as possible.

Adjusting stimulation energies


In a stimulated muscle, the number of recruited fibres depends on the stimulation energy. For programs in the Personal Trainer category, which involve powerful muscular contractions, you must therefore use maximum stimulation energies (up to 694), always at the limit that you can endure, in order to recruit the maximum number of fibres.

V HOW IT WORKS

You are strongly advised to read carefully the contraindications and safety measures described at the start of this manual in Chapter 1 prior to using your stimulator.

Description of the device

- A** On/Off button
- B** "i" button used for:
- Immediate access to information displays (illustrations of electrode placement or explanations of the programs)
 - Increasing stimulation energies in several channels simultaneously
- C** Sockets for the 4 stimulation cables
- D** Stimulation cables
- | | |
|--------------------|-------------------|
| Channel 1 = blue | Channel 2 = green |
| Channel 3 = yellow | Channel 4 = red |
- E** +/- buttons for the 4 stimulation channels
- F** Socket for the charger (Slide the red cover to the right to free the charger connector)
- G** Battery compartment
- H** Belt clip socket

EN

Battery insertion

Open the cover of the battery compartment and insert the battery with the label upwards so that the + and - poles of the battery are facing the contacts on the device. Then replace the lid. If the device is not to be used for more than 3 months, ensure the battery is fully charged. If the device is not to be used for more than 6 months, ensure the battery is fully charged and remove it from the stimulator. Before removing the battery, switch off the stimulator.

Connections

Connecting the charger

The Compex has considerable operating autonomy, as it uses a battery. To recharge, connect the charger supplied with your device to the base of the device, then plug the charger into a socket. Remove any stimulation cables connected to the stimulator before recharging it. We strongly recommend you to charge the battery fully before using it for the first time as this will improve its performance and life span.

Connecting the cables

The stimulator cables plug into the 4 sockets on the front of the device. Cables can be connected simultaneously to the four channels of the device. Both the sockets and the cables are colour-coded to simplify use and facilitate identification of the different channels:

blue = channel 1 green = channel 2 yellow = channel 3 red = channel 4

Preliminary settings

Before using the unit for the first time, you have the possibility to choose the working language of the device and to specify a certain number of settings (display contrast adjustment, backlight adjustment, and sound volume setting). This options screen can subsequently be accessed by switching off the device and then pressing the On/Off for more than 2 seconds.

- B** Use the channel 1 +/- button to choose the language you wish to use.
- C** Use the channel 2 +/- button to adjust the contrast of the screen.
- D** Use the channel 3 +/- button to adjust the volume.
- E** Use the channel 4 +/- button to adjust the backlight.
- On:** Backlight always on.
- Off:** Backlight always off.
- Auto:** Backlight activated whenever a button is pressed
- A** Press the On/Off button to confirm the selected parameters. Your stimulator saves your options. It is now ready for use with the settings you have selected.

Selecting a category

To switch on the stimulator, press the On/Off button briefly.


- A** Press the On/Off button to turn off the unit.
- B / D** Pressing +/- on channels 1 and 3 enables a category to be selected.
- E** Press the channel 4 +/- button to confirm your choice.

Personal Trainer category

The Personal Trainer category allows you to select a goal to be achieved within 6 weeks on the basis of 3 sessions per week. When choosing a goal, it is especially useful to consult chapter 6 of this manual.

On first use, you must choose a new goal.

- A** Press the On/Off button to return to the previous screen.
- B** Use the channel 1 +/- button to choose the goal of your choice.
- E** Press the channel 4 +/- button to confirm your choice.
- N.B.:** "i" button gives access to diagrams of the electrode placements and body position recommended for this goal.

Subsequently, on selecting the Personal Trainer category, you will find your goal has been saved. In this example, two goals have been preselected.

- A** Press the On/Off button to return to the previous screen.
- B** Use the channel 1 +/- button to select one of your goals.
- C** Delete the goal by pressing the +/- button of channel 2.
- D** Use the channel 3 +/- button to create a new goal.
- E** Press the channel 4 +/- button to confirm your choice.
- N.B.:** "i" button gives access to diagrams of the electrode placements and body position recommended for this goal.

When a goal has been selected, an information screen summarises the goal together with session frequency and duration as well the number of sessions accomplished.

- A** Press the On/Off button to return to the previous screen.
- E** Press the channel 4 +/- button to start the program.
- N.B.:** A session has to be fully completed in order to be recorded as having been accomplished.
- N.B.:** "i" button gives access to diagrams of the electrode placements and body position recommended for this goal.

Wellness Category


The Wellness category contains programs relating to wellbeing and pain relief. When choosing a program, it is especially useful to consult chapter 7 of this manual.

- A** Press the On/Off button to return to the previous screen.
- B** Use the channel 1 +/- button to select a program.
- E** Use the channel 4 +/- button to confirm your choices and launch the program.
- N.B.:** "i" button gives access to a short description of the program.

Adjusting stimulation energies

When you start a program, you will be prompted to increase the stimulation energies. This is critical to the success of any session. To find out what energy level needs to be attained for each program, refer to chapters 6 and 7 of this manual.


A Press the On/Off button to place the unit in Pause mode.

BCDE The Complex bleeps and the symbols of the four channels flash, changing from + to 000: the four channels are at 000 energy. You must increase the stimulation energy so that the stimulation can start. To do this, press the + buttons for the relevant channels until the desired setting is reached. If you want to increase the energy level of all four channels simultaneously, press the "I" button, located below the On/Off button. Press the "I" button twice to increase the levels in the first 3 channels, and 3 times to increase the levels in the first 2 channels. When you activate the "I" button, a black bar highlights the combined channels.

Program progression

Stimulation actually starts when the stimulation energy has been increased. The examples reproduced below are sufficient to understand the general rules.

- A** Press the On/Off button to interrupt the program momentarily. To restart it, simply press the channel + for one of the channels. The session will resume at 80 % of the energy levels that were being used prior to the interruption.
- B C D E** The different energies reached during the contraction phase are shown by a series of black bar graphs, energies of the rest phase are shown by hatched bar graphs.
- N.B.:** Active rest phase stimulation energies are automatically set at 50% of contraction energies. These can be modified during the rest phase. Once modified, they will be totally independent of the contraction energies.

EN

End of program

At the end of each session, the following screen is displayed and a short melody plays. To switch off the stimulator, press the On/Off button twice.

Battery Level and Recharging

Never recharge the stimulator without first disconnecting the stimulation cables. Always use the charger provided by Compex to recharge the batteries.

The Compex runs on a battery. The battery's performance depends on the programs and stimulation energy used.

We strongly recommend you to charge the battery fully before using it for the first time as this will improve its performance and life span.

If the device is not to be used for more than 3 months, ensure the battery is fully charged. If the device is not to be used for more than 6 months, ensure the battery is fully charged and remove it from the stimulator. Before removing the battery, switch off the stimulator.


Battery level

The battery icon on the top right of the screen indicates the level of charge in the battery. The battery icon flashes when the battery is totally flat. It is no longer possible to use the device. Recharge it immediately.

Recharging

Remove all stimulation cables from the stimulator before recharging it. Plug the charger into a socket and connect the stimulator by sliding the red cover to the right to free the charger connector.

The charge menu illustrated below appears automatically.


The duration of the charging operation is shown on the screen (a complete charge may take 10 hours). The battery icon is animated while the battery is recharging.

When fully charged, the icon will be full and the total time taken to recharge the battery will flash on the screen. Simply disconnect the charger: the Compex will turn off automatically.

Problems and solutions

Electrode fault

The Compex "bleeps" and alternatively displays the symbol of a pair of electrodes and an arrow pointing to the channel where a problem has been detected. In the example above, the stimulator has detected an error in channel 1.

- Check that electrodes are connected to this channel.
- Check whether the electrodes are old, worn, and/or the contact is poor: try using new electrodes.
- Try using the stimulation cable on a different channel. If the cable is still showing a fault, replace it (www.cefarcompexstore.com).


Stimulation is not producing the usual sensation

- Check that all the settings are correct and ensure the electrodes are positioned properly.
- Change the positioning of the electrodes slightly.

The stimulation effect causes discomfort

- The electrodes are beginning to lose adhesion and no longer provide good contact on the skin.
- The electrodes are worn and need to be replaced.
- Change the positioning of the electrodes slightly.

The stimulator is not working

- If an error screen appears while the equipment is in use, note the error number (in the example error 1/0/0) and contact the consumer service listed and approved by Chattanooga Group.

VI PERSONAL TRAINER GOALS

You wish to firm your thighs

The thigh muscles are among those placed under the greatest stress during locomotory movement such as walking, and to an even greater extent during physical activities that involve running, jumping, movement through water, etc. A sedentary lifestyle impairs the quality of these muscles, so that they become increasingly slack and no longer perform their supporting function correctly, which results in the thighs having a flabby and plump appearance.

The goal **Firm thighs** has the objective of submitting the muscles of progressively greater levels of activity, thus enabling restoration and/or maintenance of a muscle tone that can improve the thigh firmness and appearance.

Length: 6 weeks, 3 sessions/week. Having completed the goal, you may continue with 1 session/week as maintenance.

Electrode position: 1

Stimulation energy: maximum bearable energy.

You wish to reduce your waist

To maintain or restore a satisfactory body image, stimulation of the abdominal muscles is of particular relevance. In fact, the tone of the muscle fibres of which these muscles are formed normally enables the muscles to support and hold the organs in place. When the tone becomes inadequate, as is generally the case when these muscles are not stimulated sufficiently due to lack of appropriate physical activity, their supporting role is not properly assured, which has negative consequences for the figure, such as the appearance of rolls of fat around waist and a soft and sometimes protruding stomach.

The goal **Reduce waist** subjects the abdominal belt muscles to appropriate and progressive stimulation, thus enabling a slimmer waist and firmer stomach to be regained.

Length: 6 weeks, 3 sessions/week. Having completed the goal, you may continue with 1 session/week as maintenance.

Electrode position: 2

Stimulation energy: maximum bearable energy.

You wish to sculpt the buttocks

The buttocks are one of the so-called "problem areas" of the female body. This is due to the elasticity of the tissues and the specific structure of the adipose cells. These two phenomena favour the accumulation of fat and water, resulting in flabbiness and an orange-peel appearance. The buttock muscles alone occupy the whole of the posterior region of the pelvis, and it is their tone which enables flabbiness of the buttocks and the resulting negative appearance of the figure to be avoided.

The goal **Sculpt buttocks** consists of making the buttock muscles work and subjecting them to the specific activity required to improve their firmness and give the buttocks a rounded and sculpted appearance.

Length: 6 weeks, 3 sessions/week. Having completed the goal, you may continue with 1 session/week as maintenance.

Electrode position: 3

Stimulation energy: maximum bearable energy.

You wish to tone your arms

In women, the anatomical region of the arms, particularly the triceps muscles, is a locations where fat readily accumulates and the skin easily becomes saggy. This sagging is in part due to muscular insufficiency and to the particular thinness of the skin.

The natural ageing of the skin, which is very much promoted by prolonged and repeated exposure to the sun, is associated with a decrease in the control exerted by the biceps and triceps muscles resulting in a progressive sagging of the tissues and a flabby appearance that is most frequently visible on the posterior part of the arms.

The goal **Tone arms** offers a progressive stimulation of the arm muscles, which is able to restore good muscle tone.

Length: 6 weeks, 3 sessions/week. Having completed the goal, you may continue with 1 session/week as maintenance.

Electrode position: 4

Stimulation energy: maximum bearable energy.

You wish to develop your shoulder muscles

The profile of the shoulders is primarily due to the deltoid muscle, the volume of which constitutes the visible rounded part of this region. The volume of a muscle is developed by making the muscle do intensive work of sufficient quantity to cause the muscle fibres to enlarge.

The goal **Tone shoulder** enables the deltoid muscles to be made the targets of high intensity work so as to obtain a visible and harmonious enlargement of the volume of the shoulder muscles.

Length: 6 weeks, 3 sessions/week. Having completed the goal, you may continue with 1 session/week as maintenance.

Electrode position: 5

Stimulation energy: maximum bearable energy.

You wish to optimise the anti-cellulite action

Even slim people, who have a firm body, may have local accumulations of cellulite, namely areas of skin characterised by infiltration of the tissue and an "orange peel" appearance. Numerous treatments for this problem exist. However, apart from surgery, none of these treatments have shown a truly satisfactory efficacy. On the other hand, a combination of local skin treatments (ultrasound, drainage, mesotherapy, suction, compression garments, etc) appears to give better results. It is solely from this viewpoint, that the goal **Prevent cellulite** is of benefit. Used alone, it cannot get rid of your cellulite, although it is clearly more effective than old fashioned electrolipolysis.

The goal **Prevent cellulite** aims to create electrical stress and vasodilation in fat cell masses and areas of cellulite, which will reinforce the action of other local treatments that are being used in combination.

Length: 6 weeks, 3 sessions/week. Having completed the goal, you may continue with 1 session/week as maintenance.

Electrode position: 6 (hips), 7 (thighs), 8 (buttocks), 9 (stomach)

Stimulation energy: maximum bearable energy.

VII WELLNESS PROGRAMS

Relaxing massage

This program is indicated for reducing unpleasant muscle tension and creating a relaxing sensation of well-being, in particular after a tiring day. It induces a state of optimal, very pleasant relaxation.

You want to get rid of uncomfortable muscle tension in the nape of the neck

Remaining seated for long periods, associated with repetitive arm movements (as is frequently the case in front of a computer screen), can provoke increased tension in the muscles in the nape of the neck. Any other condition of stress can also lead to a state of excessive muscle tension, which is often responsible for painful or uncomfortable sensations. The in-depth effects produced by the **Relaxing massage** program ease such painful sensations, with a particularly relaxing effect.

Usage: to be used on the muscles in the nape of the neck, whenever an occasional sensation of painful muscular tension is felt; to be repeated, if necessary, in cases of particularly severe muscle tension.

Electrode position: 14

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches.

Anti-stress massage

This program enables a return to a feeling of well-being after stressful situations. It promotes the reduction of stress-related muscle tension.

Returning to a state of well-being after stressful situations

The many stresses of daily life and the worries that may accompany them cause a certain degree of mental tension to which everyone reacts differently. This hyperactivity of the brain and the emotions also has repercussions for our bodies, which may respond in many unwanted ways.

Undoubtedly, the most frequent problem is a sensation of discomfort associated with persistent muscular tension. This is sometimes intensified by physical constraints, such as the need to remain for a prolonged periods in positions unsuited to the well-being of our body. For a person under stress, it is difficult to relax. The muscles are always too tense, producing uncomfortable and even painful sensations. In this kind of situation, the **Anti-stress massage** program can reduce muscle tension by means of low-frequency stimulation, thus producing a sensation of perfect calm. Applied to the affected region, at the rate of two sessions a day, this program eliminates the unpleasant sensations produced by a stressful situation and produces a relaxing effect leading to a return to a state of well-being.

Usage: use this program daily during the stressful period or whenever a painful muscular sensation is felt.

Electrode position: 14

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches.

Undulated massage

The aim of this program is to eliminate uncomfortable physical tension and to regain or achieve a state of optimal relaxation. It uses a stimulation mode that produces a particularly pleasing undulating effect, which increases the efficacy of this program.

To eliminate unpleasant physical sensations and recover an optimal state of relaxation

To relax and eliminate uncomfortable physical sensations after a difficult and often wearisome day is the goal of many. The **Undulated massage** program offers an ideal solution for anyone who wants to improve their well-being in this way. The automatic succession of the different program sequences allows the user to obtain beneficial selective effects to overcome the various habitual components of physical discomfort: circulatory slow-down, exaggerated muscle tone, micro-stresses on the joints, etc. Also, the wave-like effect obtained through the alternation of the stimulation phases on the 4 channels further increases the efficacy of the program and produces a pleasant sensation that is particularly appreciated. The program can be used effectively on all parts of the body; however, the back and its muscles, which support the spine, are particularly prone to these problems.

Usage: use this treatment at will, whenever you feel a desire and/or need.

Electrode position: 18

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches.

Practical guidelines for use

The **Undulated massage** program uses a new stimulation mode, which produces a particularly pleasant wave-like effect that increases the efficacy of this program. In most of the sequences of this program, the stimulation is alternated among the four channels giving a faithful reproduction of various techniques of manual massage, which must always be performed in the direction of the venous return flow, that is to say from the extremities toward the heart. This is especially important in order to obtain the most benefit from the circulatory effects of these techniques and programs. To benefit fully from this wave-like effect, particular care should be given to the placement of the electrodes, and to the rules for their use as indicated below.

Rule no. 1: number of stimulation channels to use

The correct use of the Undulated massage program requires 4 stimulation channels, and therefore 8 electrodes, to be used.

These 4 stimulation channels must be positioned:

- on a limb or part of a limb (e.g.: a thigh),
- on the back,
- on two limbs (e.g. : both calves).

Rule no. 2: required order of the stimulation channels

The wave-like effect must follow the direction of the venous return flow. That is why the required order of the four stimulation channels must be strictly complied with. The order of the stimulation channels is as follows::


Channel 1: cable fitted with a blue colour indicator,

Channel 2: cable fitted with a green colour indicator,


Channel 3: cable fitted with a yellow colour indicator,

Channel 4: cable fitted with a red colour indicator.


It is always recommended to follow the general indications below:

Application on a limb or part of a limb

Application on the back

Application on both limbs

Heavy legs

This program is intended to eliminate the sensation of heaviness and discomfort in the legs which is felt in particular circumstances (standing for prolonged periods, excessively high ambient temperature, hormonal imbalance related to the menstrual cycle, etc.). It quickly induces a feeling of lightness and suppleness.

You want to prevent the feeling of heaviness and discomfort in your legs

The feeling of heaviness in the legs is due to a temporary deficiency in the venous return. The insufficient supply of oxygen to the tissues, particularly to the muscles, that results from this situation, leads to this uncomfortable sensation in the legs. Use the **Heavy legs** program to accelerate the venous return and obtain a significant relaxing effect on painful muscles.

Usage: to be used whenever an occasional feeling of heaviness in the legs occurs.

Electrode position: 16

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches.

Lymphatic drainage

The deep massage produced by this program boosts the circulation of blood and fluids through the tissues. It thus enables the elimination of excess water and accumulated toxins, as well as regeneration and cleansing of the skin tissue.

Treatment of occasional swelling of the feet and ankles

In unusual circumstances, impaired circulation of the venous blood often leads to an accumulation of blood and lymph in the lower limbs. This phenomenon, which manifests as swelling of the ankles and feet, causes a disagreeable sensation of heaviness and tension in the affected areas. The deep massage produced by the **Lymphatic drainage** program is effective in boosting the return lymphatic circulation and thus promoting the resolution of these problems.

Usage: for 8 weeks, 3 sessions/week. Use this programme when there is occasional swelling of the feet and ankles.

Electrode position: 20

Stimulation energy: increase the stimulation energy progressively until very pronounced muscle twitches occur during the active rest phase and powerful tetanic contractions during the contraction phase.

Pain management TENS

All painful phenomena, especially pain in and around the joints, can be effectively relieved using this program, which blocks the transmission of pain by the nervous system.

You have pain in the region of the knee

The knee is one of the joints that is most subjected to stress in the course of physical and sports activities, and also during simple movements in daily life, such as walking and going upstairs, for example. Wear of the cartilage, stretching of a ligament or irritation of a tendon may cause knee pain, which can be effectively relieved by rest or by use of the **Pain management TENS** program.

Usage: use in accordance with the course of the pain; several sessions per day may be performed if required.

Electrode position: 21

Stimulation energy: gradually increase the stimulation energy until a pronounced tickling sensation can be felt under the electrodes.

Caution: Do not use this program for prolonged periods without medical advice.

Chronic pain

Thanks to its analgesic action, this program provides rapid relief in case of chronic muscle pain (back pain, lumbar pain, etc.). It boosts the circulation, oxygenates and relaxes stiff muscles.

You have muscular tension in the nape of the neck, middle back or lower back region

Pain in the back or nape of the neck is due to chronic contraction of the muscles, with accumulation of toxins. The pain can then become continuous or appear after only a few minutes spent working in an unfavourable position (for example, working in front of a computer screen). This chronic pains can be treated effectively with the **Chronic pain** program which, in addition to relieving pain, reactivates the circulation and relaxes the contracted muscles.

Usage: for 4 weeks, 2 sessions/day, with a 10-minute break between the 2 sessions.

Electrode position: 14 (nape of neck), 13 (middle back) or 12 (lower back)

Stimulation energy: increase the stimulation energy gradually until you obtain clearly visible muscular twitches.

Caution: Do not use this program for prolonged periods without medical advice.

VIII EMC TABLES

The Compex Fit needs special EMC precautions and must be installed and started according to the EMC information supplied in this manual.

All RF wireless transmission systems can affect the Compex Fit.

The use of accessories, sensors and cables other than those recommended by the manufacturer may result in stronger emissions or reduce the immunity of the Compex Fit.

EN

The Compex Fit should not be used beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the Compex Fit works properly in the chosen configuration.

RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER CONCERNING ELECTROMAGNETIC EMISSIONS

The Compex Fit is intended for use in the electromagnetic environment specified below.

The customer or user of the Compex Fit should ensure that it is used in this environment.

Emissions test	Compliance	Electromagnetic environment - guide
CISPR 11 RF emissions	Group 1	The Compex Fit uses RF energy only for its internal operation. Consequently, its RF emissions are very low and are unlikely to interfere with any adjacent electrical device.
CISPR 11 RF emissions	Class B	The Compex Fit is suitable for use in any establishment, including a private dwelling and a place connected directly to the low voltage mains supply which powers residential buildings.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ emission oscillations IEC 61000-3-3	Not applicable	

RECOMMENDED SPACING BETWEEN A PORTABLE AND MOBILE COMMUNICATION APPLIANCE AND THE COMPEX FIT

The Complex Fit is designed for use in an electromagnetic environment in which radiated RF waves are controlled. The buyer or user of the Complex Fit can contribute to preventing electromagnetic interference by maintaining a minimum distance between RF portable and mobile communication appliances (transmitters) and the et Complex Fit according to the table of recommendations below and according to the maximum output power of the telecommunication appliance.

Maximum transmitter output power W	Spacing according to the frequency of the transmitter m CISPR 11		
	From 150 kHz to 80 MHz	From 80 MHz to 800 MHz	From 800 MHz to 2.5 GHz
	$d = 1.2 \cdot \sqrt{P}$	$d = 1.2 \cdot \sqrt{P}$	$d = 2.3 \cdot \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

In the case of transmitters whose maximum output power is not shown in the table above, the recommended spacing of d metres (m) can be calculated using the appropriate equation for the transmitter frequency, where P is the maximum output power of the transmitter in watts (W) as set by the transmitter manufacturer.

NOTE 1: At 80 MHz and at 800 MHz, the spacing for high frequency amplitude is applied.

NOTE 2: These guidelines may not be appropriate for some situations. Electromagnetic wave propagation is modified by absorption and reflection due to buildings, objects and persons.

RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER - ELECTROMAGNETIC IMMUNITY

Compex Fit is designed for use in the electromagnetic environment stipulated below. The buyer or user of the Compex Fit must ensure it is used in this recommended environment.

Immunity test	Test level IEC 60601	Observance level	Electromagnetic environment - guide
Electrostatic discharge (DES) CEI 61000-4-2	± 6 kV at the contact ± 8 kV in air	± 6 kV at the contact ± 8 kV in air	Floors must be wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity must be maintained at a minimum of 30%.
Rapid electrical transients in bursts of noise CEI 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable System battery-powered	The quality of the power supply should be that of a typical commercial or hospital environment.
Surge currents CEI 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Not applicable System battery-powered	The quality of the power supply should be that of a typical commercial or hospital environment.
Voltage trough, short power cuts and voltage variations on power supply input lines CEI 61000-4-11	< 5% U_T (trough >95% of U_T) for 0,5 cycle 40% U_T (trough 60% of U_T) for 5 cycles 70% U_T (trough 30% of U_T) for 25 cycles <5% U_T (trough >95% of U_T) for 5 sec	Not applicable System battery-powered	The quality of the power supply should be that of a typical commercial or hospital environment. If the Compex Fit user requires continuous operation during mains power cuts, we recommend that the Compex Fit is powered by a UPS or a battery.
Magnetic field at mains frequency (50/60 Hz) CEI 61000-4-8	3 A/m		Magnetic fields at the mains frequency should be at the level of a representative site located in a typical commercial or hospital environment.

NOTE: U_T is the AC supply voltage before application of the test level.

Immunity test	Test level IEC 60601	Observance level	Electromagnetic environment - recommendations
			Portable and mobile RF communication devices must only be used relative to the Complex Fit and its wiring at a distance which is not less than the spacing recommended and calculated using the appropriate equation for the transmitter's frequency. Recommended spacing
RF immunity IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Signal lines Not applicable for power supplies of battery-powered systems	$d = 1,2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz 10 V/m 26 MHz to 1 GHz	3 A/m 10 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power of the transmitter in watts (W) set by the manufacturer's specifications and where d is the recommended spacing in metres (m). The field intensity of RF fixed transmitters, as determined by an electromagnetic survey a must be less than the observance level to be found in each frequency range b. Interference may occur close to any appliance identified by the following symbol: 

NOTE 1: From 80 MHz to 800 MHz, the spacing for high frequency amplitude is applied.

NOTE 2: These guidelines may not be appropriate for some situations. Electromagnetic wave propagation is modified by absorption and reflection due to buildings, objects and persons.

^a The field intensity from fixed transmitters, such as radio telephone base stations (cellular/wireless) and a mobile radio, amateur radios, AM and FM radio emissions and TV emissions cannot be predicted with any accuracy. It may therefore be necessary to consider an analysis of the electromagnetic environment of the site to calculate the electromagnetic environment coming from RF fixed transmitters. If the field intensity measured in the environment where the Complex Fit is located exceeds the appropriate RF observance level above, the Complex Fit should be monitored to ensure it is operating properly. In the event of abnormal operation, new measures may then be imposed, such as realignment or movement of the Complex Fit.

^b Above the frequency amplitude from 150 kHz to 80 MHz, the field intensity must be 3 V/m