**MANUALE D'USO**

**USER MANUAL**

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**MANUAL DE USO**

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**ELECTROSTIMULATOR**

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# ENGLISH

- Contents	29
- Equipment	30
- Contraindications – Warnings – Maintenance	30
- Technical features	33
- Method of use	34
- Description of the device	34
- Stimulation intensity	35
- Description of treatments	36
- Programmability	39
- Ionophoresis	40
- ELITE 2 PROGRAMME LIST	42
- GENESY 2 PROGRAMME LIST	50
- REHAB 2 PROGRAMME LIST	53
- Frequently Asked Questions	55
- Guarantee	56
- Purchasable accessories	56
- Electrode Placement	141

### Equipment:

- a. 2 coloured cables for electrode connection
- b. Bag containing 4 reusable self-adhesive square electrodes
- b. Bag containing 4 reusable self-adhesive rectangular electrodes
- d. Carrying bag
- e. User manual with electrode placement guide
- f. Unit
- g. Guarantee slip
- h. CD ROM Personal Trainer only in mod. ELITE 2
- i. 4 non-rechargeable AA batteries

### Purchasable accessories

- Probes for urology treatments
- Special electrodes for ionophoresis

The electrostimulation device is supplied with cables and electrodes for use. Once the package has been opened, check that the basic equipment is complete. If any element is missing, immediately contact the authorized dealer from which you purchased the product.

Visually check the integrity of the device and electrodes.

**Electrodes:** The supplied electrodes are for use on a single patient, self-adhesive and pregelled with a female cable. They are EC-marked according to directive 93/42/EEC concerning medical devices.

All the information provided may be subject to change without advance notice.

### STANDARD SYMBOLS

The marks indicate that, according to EC regulations, the product is a class II and type BF medical device .

Caution: read the instructions contained in the attachment

This symbol on your device indicates that it meets the requirements of directive (93/42/EEC) concerning medical devices. The number of the notified body is 0476.

BF applied parts

Insulation class II

On/Off

### CONTRAINDICATIONS AND WARNINGS

Should not be used in the following cases:

- Stimulation of the front part of the neck (carotid sinus).
- Pacemaker patients.
- Tumour patients (see your oncologist).
- Stimulation of the brain region.
- Pain of unknown aetiology.
- Skin sores and other dermatological pathologies.
- Acute trauma.
- Stimulation on recent scars.
- Pregnancy.
- It is absolutely forbidden to use the electrostimulator in the eye or genital region.
- Use of the electrostimulator in water (the device is not water-proof) and in atmospheres with combustible gases, flammable vapours or in oxygen chambers.
- Use of the device on or near areas characterised by metallic synthesis means and intratissular metals (such as prostheses, osteosynthesis materials, spirals, screws, plates), if single-phase currents are used such as interferentials and direct current (ionophoresis);
- Caution is advised in case of patients suffering from capillary fragility, as excessive stimulation may cause capillary rupture.

### **Possible undesired effects**

Isolated cases of skin irritation may occur in subjects with particular skin sensitivity. In case of allergic reaction to electrode gel, suspend the treatment and consult a specialist.

Should symptoms of tachycardia or extrasystole occur during the treatment, suspend the electrostimulation and seek medical advice.

### **Particular precautions for Urology treatments**

- Patients with extra-urethral incontinence must not be treated with the stimulator.

- Patients suffering from excessive incontinence due to passage problems must not be treated with the electrostimulator.

- Patients with severe urinary retention in the upper urinary tract must not be treated with the stimulator.

- Patients with complete peripheral denervation of the pelvic floor must not be treated with the stimulator.

- Patients suffering from a total/subtotal prolapse of the uterus/vagina must be stimulated with extreme caution.

- Patients with urinary tract infections must be treated for these symptoms before stimulation.

- Before removing or touching the probe it is necessary to turn the stimulator off or adjust the intensity of both channels to 0.0 mA.

- The treatment is a personal medical prescription: do not lend the stimulator to other people.

### **Maintenance and cleaning of the device**

- In case of deterioration of the device casing, the latter should be replaced and no longer used.

- In the event of a real or presumed breakdown, do not open or to try to repair the device. Do not operate on the machine; do not open it. Repairs can be carried out only by an authorised specialised centre.

- Avoid violent impacts that could damage the device and cause malfunctions, possibly not immediately detectable.

- Keep the electrostimulator far from water sprays, heat sources, dust, damp environments and direct exposure to solar rays and high electromagnetic emission sources. Do not wrap the device in cloth or anything else during use.

- Do not clean the electrostimulator with diluents, gas oil, kerosene, alcohol or chemical products. For the cleaning, a slightly damp soft cloth is sufficient for both the device and the cables.

- After a long period of non-use, check the proper functioning of the device and batteries (apply to an Authorised Service Centre).

- The device must not be used when being charged.

- The device must be connected to the electric line through its own power supply. Before performing this operation, check that the system is in conformity with existing national regulations.

- The use of the device simultaneously with other electronic devices is inadvisable, particularly with those that are used to maintain vital functions.

### **Conservation and use of electrodes and cables**

If the cables or electrodes are deteriorated, they should be replaced and no longer used.

Clean the skin carefully before applying the electrodes.

After using the single-patient multi-use and/or mono-use electrodes, keep them with their own plastic film in a clean, closed plastic bag.

Prevent the electrodes from touching or overlapping.

Once the package has been opened, the electrodes can be used for 25-30 applications.

The electrodes should, however, always be replaced if they are not perfectly adherent to the skin.

If using non-self-adhesive electrodes, clean their surface with water.

Make sure the cables and electrodes are clean before every application.

In case of deterioration of the cables, these should be replaced with original spare parts and no longer used.

Use only electrodes supplied by the device manufacturer.

### **Warnings**

- The device must be used exclusively by mentally competent adults.

- The electrostimulator must not be used for purposes other than transcutaneous neurostimulation.

- Neurostimulation devices must be used following the instructions and under the supervision of a doctor or qualified physiotherapist.

- The electrostimulator must be used with suitable transcutaneous neurostimulation electrodes.

- The electrostimulator must be kept out of children's reach.

- The electrostimulator's current can cause interference with ECG monitoring devices.

- Do not apply current in trans-thoracic mode to prevent arrhythmia and interference with heart frequencies. (Do not perform the treatment on pectoral and dorsal muscles simultaneously).

- Do not use the device near short-wave therapy devices.

- In case of health problems, seek medical advice before using the device.

### **Intended use**

Carefully read the user manual to avoid incorrect installation, use and maintenance operations which may damage the device, the operator or the user.

Electrostimulators are designed to be used in the following environments:

- home;

- clinical;

- physiotherapy;

- rehabilitation;

- pain treatment;

- any type of sequential electrostimulation;

- cosmetic and sports purposes.

The use of this device is permitted to patients (provided that they are adequately informed of the techniques used by the device), medical and paramedical personnel.

### **Disposal of the device**

Do not throw the device or parts of it into the fire, but dispose of the product in specialised centres and in compliance with current national regulations.

## TECHNICAL FEATURES

### *EMS and TENS:*

Outputs available: Channels 1-2  
Const. current: Yes  
Intensity: 0-100 mA with  
1000 Ohm load  
Waveform: rectangular, biphasic,  
symmetric, compensated  
Frequency range: 1-120 Hz  
Pulse width: 50-400 microseconds  
Work time: 1 to 60 seconds  
Recovery time: 0 to 1 minutes  
Power supply: 4 1.5 V AA LR6 batteries

### *Ionophoresis (Genesy 2 PRO, Rehab 2 PRO models):*

Outputs available: Channel 1  
Const. current: Yes  
Minimum intensity: 0 mA/1000 Ohm  
Maximum intensity: 5 mA/1000 Ohm  
Minimum time: 1 minute  
Maximum time: 99 minutes

### *Storage and transport conditions:*

Temperature: 10-40°C  
Air humidity: 20-65%  
Atm. pressure: 700 hPa - 1060 hPa

## **Battery installation**

If necessary, switch off the device by pressing the ON/OFF button for 2 seconds. Remove the battery cover located on the back of the unit by twisting off the screw, and access the battery compartment. Insert 4 x 1.5V AA LR6 Alkaline batteries (or AA rechargeable batteries), replace the cover and twist on the screw. The batteries must be positioned with the correct polarity indicated on the compartment. An incorrect positioning of the batteries may damage the unit.

## **Cables connection**

It is possible to use one or two cables depending on how many channels (one or two) are used. For the connection, insert the cable connector into the appropriate lodging located on the upper part of the unit. The lodgings are located

precisely under the corresponding channel (left or right). Insert the cable under the selected channel. **For IONOPHORESIS use Channel 1 only.**

## **Application of electrodes**

Take the electrodes out of the original package; all the new electrodes have a sealed package. Make sure that the device is turned off.

Begin by connecting the two cable jacks to the electrodes, then remove the electrodes from their seat and place them onto the skin. For the correct positioning of electrodes, refer to the pictures provided by this manual.

After use, replace the electrodes in the appropriate seat.

**CAUTION:** Do not remove the electrodes if the unit is working.

## **WARNINGS AND INFORMATION FOR FIRST USE**

Before using the device for the first time, it is necessary to:

- Read this user manual carefully;
- Visually examine the device and all its parts to identify any signs of external damage;
- If all the segments appear on the display when the device is turned on, turn it off and then again on; if all the segments appear again, contact the customer service and do not use the device.
- If the software version/model appear on the display when the device is turned on, the latter is ready for use.
- If the machine cuts out immediately after turning it on, this means that the battery level is insufficient for use. In this case, replace the batteries following the instructions in the "Battery installation" section.

## METHOD OF USE


In order to connect the product to the user correctly, proceed as follows:

- connect the cables to the device cable sockets;
- connect the electrodes to the relative cables connected to the device;
- apply the electrodes to the skin.

The device and its accessories must be replaced in the appropriate bag in compliance with the following storage and transport conditions.

Temperature: 5-45 °C  
Air humidity: 20-85 %

## DESCRIPTION OF THE DEVICE

- A - Channel 1 intensity adjustment buttons
- B - Channel 2 intensity adjustment buttons
- C - Phase skip buttons (if present)
- D - Buttons for programme number choice, time modification and selection of parameters
- E - ON/OFF - OK - PAUSE button
- F - Function Key (Fn)

## CAUTION !

*For the ionophoresis treatment, use Channel 1 with the specific electrodes only. Electrodes are available upon payment.*

### Turn-on

To turn on the device, press the ON/OFF (OK) button for three seconds until the acoustic signal is heard. The software version/model is displayed for three seconds with a number in the lower left display.

### Programme number choice.

Select the programme number using the P+ and P- buttons. Press Ok to confirm.

During the programme number choice, the display shows the number of phases in the programme (top) and the programme total duration (top right).

### Start of the programme.

Once the programme has been selected, increase channel intensity by using the channel intensity adjustment buttons.

### Display during execution.

During execution, the display shows the programme number (top left), the total number of phases and the phase in progress (full phase) and the remaining time of the phase in progress (top right).

### Programme pause.

To pause the programme, press the OK button on the Joypad. The time stops. Press OK again to resume to the programme. Resumption after a pause repeats the last interrupted cycle from the beginning with set current values. During pause, the blinking wording PAUSE appears on the display.

### Programme interruption.

Should it be necessary to interrupt the programme, turn off the device by pressing the OK button for three seconds.

### **Time display.**

During programme execution, the display shows the remaining time of the phase in progress.

### **Programme Time, Frequency and Width Modification.**

Press the Function Fn key several times to select the parameters.

#### **Pressing the Fn key 1 time: phase duration modification.**

Once the programme is started, the phase duration can still be modified. To change the phase duration, press the Fn key; the time of the phase in progress begins to blink.

Modify the time by pressing the P+ or P- buttons. The set time will automatically be confirmed after being displayed for 5 sec.

#### **Pressing the Fn key 2 times: frequency modification of the programme in progress.**

The second time the Fn key is pressed, the word **FREQ** is displayed and the value starts to blink.

Change the frequency by pressing the P+ or P- buttons of the Joypad. The set frequency will automatically be confirmed after being displayed for 5 sec.

#### **Pressing the Fn key 3 times: width modification of the programme in progress.**

The third time the Fn key is pressed, the word **WIDTH** is displayed and the value starts to blink.

Change width by pressing the P+ or P- buttons of the Joypad. The set width will automatically be confirmed after being displayed for 5 sec.

#### **Phase skipping.**

In order to skip to the following phase before the end of the phase in progress, press the button --> for three seconds. To return to the previous phase, press the button <-- for three seconds.

#### **Increase/decrease intensity.**

To increase/decrease channel intensity, press the Up or Down buttons of the corresponding channels.

### **Keypad lock.**

Once the desired intensity level has been reached, the keypad may be locked to prevent accidental use.

Simultaneously press the *Fn button and the Ch 1 intensity increase button* for 3 seconds. the lock symbol (a key) appears on the display and all the buttons are locked except the ON/OFF button. When the keys are locked, it is impossible to increase or decrease the intensity or change the programme. To unlock the keys, press the same combination of buttons (*Fn button + Ch 1 intensity increase button*) again for 3 seconds. The lock symbol (a key) disappears and all the buttons can be used, again.

### **STIMULATION INTENSITY**

#### **Open circuit**

The device is provided with a current supply control system. If the intensity value is increased to over 10 mA and the current circuit is open (cables not connected to the device and electrodes not in contact with the skin), the electrostimulator immediately brings the intensity value back to 0 mA.

Before beginning any treatment, make sure that the cables are connected to the device, that the electrodes are applied to the area to be treated and that they are in good condition so as not to prejudice their conduction capacity.


#### **Stimulation intensity**

The value of the necessary current intensity to obtain a certain contraction is determined on an individual basis, as it can depend on the positioning of electrodes, the adipose layer, sweating, hairy areas, etc.

For these reasons the same current intensity can produce sensations that vary from person to person, day to day, and left to right side. Even during the same work session, due to adaptation phenomena, it may be necessary to adjust the intensity to obtain the same contraction.

The current intensity to be used in the various phases is suggested with an indicative value, which must be taken into consideration according to individual sensations.

- Low: the muscle is not fatigued even in case of lengthy treatments; the contraction is bearable and pleasant. First level in the intensity chart.
- Intermediate: the muscle contracts visibly but causes no articular movement. Second level in the intensity chart.
- High: the muscle contracts considerably. The muscular contraction may cause the extension or flexion of the limb if the latter is not blocked. Third level in the intensity chart.
- Maximum: maximum muscle contraction. Very taxing work. Perform only after several applications.

Recommended intensity values are indicated in the treatment descriptions.

**NOTE** The current values shown are merely indicative.

## TREATMENT DESCRIPTIONS

Following is a short description of the programmes provided by the electrostimulator.

### *SPORT*

**WARM-UP:** Increases muscular temperature for better training preparation. Increases the blood flow to and metabolism of the muscle. Recommended before training sessions or competitions.

Intensity: low or intermediate.

**MAXIMAL FORCE:** Improves maximal force capability and increases muscular mass. Allows training with maximum load and reduced trauma risk as compared to classical training. Recommended for activities requiring great physical strength. Also useful during the early preparation period as it trains the muscles to receive significant stress.

Intensity: maximum, at pain threshold limit. During active Rest, adjust the intensity to obtain comfortable contractions.

**RESISTANCE FORCE:** Improves the ability to display a high level of force for an extended period of time and to resist toxin accumulation, thus delaying muscular fatigue. The programme involves a long series of lengthy contractions, alternated with brief active rest. Recommended for sports requiring intense muscular work for an extended period of time.

Intensity: intermediate or maximum. Active Rest should induce muscular vibrations.

**EXPLOSIVE FORCE:** The Explosive Force programme increases the ability to rapidly achieve maximum levels of force. The objective is to recruit the maximum number of muscle fibres in the shortest time possible. It is recommended to begin with a Maximal Force programme to prepare the muscle for the Explosive Force programme. Recommended in all sports requiring maximal force under natural load conditions,

such as sprints, ball sports, and alpine skiing.  
Intensity: maximum muscular contraction at pain threshold limit. Active Rest should induce muscular vibrations.

**RELAXATION:** Enables effective muscle relaxation thanks to the pleasant contractions that increase blood circulation. These programmes can be used whenever muscular relaxation is needed.  
Intensity: low to intermediate.

**CAPILLARISATION:** Increases the blood flow to the muscle and improves muscle resistance and recovery. An increased arterial capacity keeps the primary and secondary capillary network active, thus improving the tissue oxygenation system. This delays fatigue during very intense work. Capillarisation programmes are particularly recommended for sports requiring resistance force and endurance.  
Intensity: intermediate.

**POST-COMPETITION / TRAINING RECOVERY:** Improves recovery from muscular fatigue after training or competitions and enables easier training resumption. The stimulation varies from muscular vibrations to light contractions. This allows greater circulation and reduced muscular hardening. It also alleviates pain thanks to the release of endorphins. This programme is valid for all sports as it improves and speeds up the recovery of maximum muscular capacity. This programme is recommended 8-24 hours after training or competitions.  
Intensity: low to intermediate.

## ***FITNESS AND BEAUTY***

**DRAINAGE:** The coordinated stimulation obtained by using all channels of the electrostimulator produces an in-depth muscular massage effect. The massage activates venous and lymphatic systems by means of a mechanical effect. The purpose is to improve blood circulation and enable the toxin elimination. Intensity: low to intermediate.

**LIPOLYSIS:** The low frequency of this programme positively affects the metabolism and increases the elimination of toxins and excess fluid. This is a method for reshaping critical areas, modifying the metabolism of lipocytes. Improves blood circulation. Helps to eliminate the orange peel skin.  
Intensity: low to intermediate.

**VASCULARISATION:** Enables effective muscle relaxation thanks to the pleasant contractions that increase blood circulation. These programmes can be used whenever muscular relaxation is needed. Intensity: intermediate.

**FIRMING:** This programme is suitable for sedentary individuals with flaccid muscles that require reactivation of muscular tone. It prepares the muscles for more intense stimulations. It can be performed every day on large muscle groups such as thighs, abdominal and dorsal muscles.  
Intensity: intermediate to high.

**TONING:** The toning programme allows to develop and maintain muscular tone. This type of treatment may be combined with aerobic activity and is suitable for users of any level. The programme is especially directed at all those who, mindful of their own physical condition, wish to integrate their own regular training. The programme can be used up to three times a week.  
Intensity: intermediate to high.

**SCULPTING:** The sculpting programme is a merely aesthetic treatment and is ideal for those who wish to better define their already toned muscle system. The work performed is very intense, so slight muscle soreness may be experienced the day after the treatment. However, this effect may be avoided by performing stretching exercises before and after the session or a winding-down programme. The programme is suitable for those who wish to achieve muscular definition in those areas where traditional muscular development exercises are insufficient. Intensity: intermediate to high.

**MASS DEVELOPMENT:** A specific programme for those who wish to increase muscle mass and

volume. Intensity: maximum muscular contraction at pain threshold limit.

A muscle warm-up programme should be carried out before performing any of the following programmes: Firming, Toning, Sculpting and Mass Development.

### **HEALTH**

Uses mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. This analgesic effect occurs during programme execution and can last for a few hours. However, its duration may vary. The TENS programmes also improves circulation. Place the electrodes on the painful area. These programmes can be used in any situation requiring pain relief. Caution! TENS programmes alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice.

**TENS (Analgesic-Endorphinic):** Allows pain relief and effective muscle relaxation. TENS programmes use mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. Muscular relaxation improves blood circulation. It can be used in any situation requiring both pain relief and muscle relaxation. Caution! TENS programmes alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice.

Intensity: ANALGETIC TENS: tingling but not painful sensation. Should not induce muscular contractions. ENDORPHINIC TENS: visible muscular vibrations.

**TENS (Cervical pains, Backache, Periarthritis, Epicondylitis, Sciatic neuralgia):** Uses mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. This programme combines stimulation with light contractions that stimulate the release of endor-

phins. Caution! TENS programmes alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice. Intensity: tingling but not painful sensation.

**HEMATOMA:** Increases the absorption of the muscular hematomas due, for example, to contusions. Improves blood circulation in the area. This programme should not be used earlier than 72 hours after the muscular trauma. Always seek medical advice if further instructions are required. Intensity: intermediate to high.

### **REHAB**

**ATROPHY RECOVERY:** Resumption of physical activity after a period of interruption. After a trauma or other interruptions of physical activity, muscles may lose their tone and mass. In these cases, normal muscular capacity must be regained as soon as possible. The Atrophy Recovery programme enables training of specific muscles concerned.

It can be used in any situation requiring activity resumption. Intensity: low to intermediate.

**RE-MUSCULATION:** Increases the diameter and capacity of muscle fibres.

This programme is recommended after using the Atrophy programme for a period, as soon as the muscle shows a slight recovery of volume and tone. Intensity: Intermediate.

**MOBILIZATION:** The stimulation is aimed at neuromuscular facilitation, which especially improves the mobilization of the injured articulation. Prepares the limb to resume its normal functionality. It is recommended to use this programme after the recovery of muscular trophism. Intensity: intermediate.

**HEMIPLEGIA - PARAPLEGIA:** If the patient suffers from a perception problem with decreased attention (semi-inattention), integrating the treatment with training of this ability may be extremely useful. The treatment must not be merely passive, but the patient can take on an active role in the training. The programme can

be used for the facilitation and rehabilitation of motor faculties as well as a method to reduce spasticity.

Intensity: depending on the muscle district to be stimulated.

**UROLOGY:** Programmes recommended for the stimulation of perineal muscles.

These programmes are recommended in case of urological pain or incontinence.

Programmes to be used with probes in compliance with the law.

**NOTE:** Because it includes clinical stimulation programmes, the apparatus is classified as a medical device. For this reason, it is certified according to directive 93/42/EEC on medical devices by CERMET, Notified Body No. 0476. The certification encompasses clinical applications.

## **PROGRAMMABILITY (PRO models only)**

To access programmability after turning on the device, press the Fn key for 3 seconds. The first of the 10 free programmes available appears. Use the P+ or P- buttons to select the number of the programme you wish to create and press OK to confirm. By using the P+ and P- buttons, select the type of impulse: continuous or intermittent --- and press OK to confirm.

### **CONTINUOUS IMPULSE STIMULATION ---**

a) Use the P+ or P- buttons to set the phase time duration and press OK to confirm. **(Pt)**

b) Use the P+ or P- buttons to set the frequency range and press OK to confirm. **(Fw)**

c) Use the P+ or P- buttons to set the work impulse width and press OK to confirm. **(Aw)**

To programme only phase 1, press the Fn key for 3 seconds. The programme will thus be saved. Press Ok to begin the treatment. To programme other phases, repeat stages a), b) and c). Once all the phases have been programmed, press the Fn key for 3 seconds.

### **INTERMITTENT IMPULSE STIMULATION**

d) Use the P+ or P- buttons to set the phase time duration and press OK to confirm. **(Pt)**

e) Use the P+ or P- buttons to set the frequency range and press OK to confirm. **(Fw)**

f) Use the P+ or P- buttons to set the rest frequency and press OK to confirm. **(Fr)**

g) Use the P+ or P- buttons to set the work/rest impulse width and press OK to confirm. **(Aw)**

h) Use the P+ or P- buttons to set the work time and press OK to confirm. **(wt)**

i) Use the P+ or P- buttons to set the rest time and press OK to confirm. **(rt)**

To programme only phase 1, press the Fn key for 3 seconds. The programme will thus be saved. Press Ok to begin the treatment. To programme other phases, repeat stages d), e), f), g), h) and i). Once all the phases have been programmed, press the Fn key for 3 seconds.

*For the ionophoresis treatment use Channel 1 with the specific electrodes only. Electrodes are available upon payment.*

*Given the medical use of these applications, they should only be performed upon consent of medical personnel.*

## IONOPHORESIS

The galvanic current maintains a constant intensity over time and allows the movement of charged particles through the tissues. If charged particles are drugs, then the direct current acts as a carrier enabling the penetration of medicinal substances into the organism. This technique is called IONOPHORESIS.

Ionophoresis can be used to administer water-soluble medicines through the epidermis. Ionophoresis may be used in all treatments that positively affect local inflammatory states such as bursitis, tendinitis, fascitis, myositis and capsulitis.

### Precautions

If the procedure followed is not rigorous and the intensity is too high, this may cause skin irritation or burns. The recommended procedure is 0.2 mA per cm<sup>2</sup> of electrode. The skin must not have any type of lesions (do not shave before treatment in order to avoid possible microcuts caused by the razorblade).

The use of IONOPHORESIS programmes is prohibited in subjects with metallic implants. Make sure that the patient is not in contact with metallic objects (e.g., table or chair).

### Use of IONOPHORESIS

The duration of the IONOPHORESIS programme is set by default to 20 minutes.

The intensity value is adjustable from a minimum of 1 mA to a maximum of 5 mA.

The safety rule for the intensity is 0.2 mA per cm<sup>2</sup> of electrode.

## Drugs for IONOPHORESIS

*Before use, always seek medical advice.*

Carefully read the instructions for use on the drug packaging.

Before any therapy, check drugs suitability by consulting a doctor or a physiotherapist, and identify their exact polarity.

The drug used for the therapy **SHOULD NEVER BE APPLIED DIRECTLY ONTO THE SKIN**, but always onto the absorbent surface of the electrode corresponding to the drug polarity, while the absorbent surface of the other electrode should be dampened with slightly salted water, to favour current circulation.

$$\text{DENSITY} = \frac{I \text{ (Displayed current)}}{\text{Electrode area}}$$

For safety **DENSITY** < 0.2 mA/cm<sup>2</sup>

## DRUGS FOR IONOPHORESIS

Some drugs for **IONOPHORESIS** are provided below. Check the correct drug polarity. Drug polarity is indicated on the drug packaging.

DRUG	ACTION AND USE
Zinc chloride (+)	Infection of maxillary sinuses, disinfection of the body
Carbaine (+)	Local anesthetic analgetic
Silver nitrate (+)	Anti-inflammatory <b>NOTE: At low current intensity.</b>
Aconitine nitrate (+)	Analgetic, trigeminal neuralgia
Hyaluronidase (+)	Anti-edematous, thrombophlebitis
Benzidamine (+)	Contusions, sprains, post-fracture after-effects, mioarthralgia, osteoarthritis
Thiomucase (-)	Backache, cervical pains, cellulitis
Sodium chloride (-)	Sclerotic and softening <b>NOTE: At low current intensity.</b>
Calcium chloride (+)	Recalcifying <b>NOTE: Do not use in case of arteriosclerosis</b>
Alphachymotrypsin (+)	Anti-inflammatory, contusions, sprains, traumatic oedemas
Lithium salicylate (+)	Analgetic, articular rheumatism
Histamine dichloride (+)	Rheumatic pain, muscular spasms, rheumatic myalgia
Procaine (+)	Contusions, analgetic, sprains, post-fracture after-effects
Potassium citrate (+)	Anti-inflammatory
Sodium salicylate (-)	Acute articular rheumatism, myalgia
Potassium iodide (-)	Arthritis, arthrosis

# PROGRAMME LIST ELITE 2 - SPORT MAN

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <b>Intensity</b>
Sport M. <b>1</b>	Capillarisation	As shown in the pictures at the end of this manual	20' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Sport M. <b>2</b>	Warm-up	As shown in the pictures at the end of this manual	10' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Sport M. <b>3</b>	Recovery, post-competition/training	As shown in the pictures at the end of this manual	20' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Sport M. <b>4</b>	Maximal force, lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>5</b>	Resistance force, lower limbs	I1-I2-I3-I4-I5-I6	5' 20' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>6</b>	Explosive force, lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>7</b>	Reactivity, lower limbs	I1-I2-I3-I4-I5-I6	5' 20' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>8</b>	Relaxation, lower limbs	I1-I2-I3-I4-I5-I6	20' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Sport M. <b>9</b>	Maximal force, abdominals and dorsals	T1-T2-T3-T4-T5-T6	5' 15' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>10</b>	Resistance force, abdominals and dorsals	T1-T2-T3-T4-T5-T6	5' 20' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>11</b>	Relaxation, abdominals and dorsals	T1-T2-T3-T4-T5-T6	20' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Sport M. <b>12</b>	Maximal force, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>13</b>	Resistance force, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 20' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>14</b>	Explosive force, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' <b>Int Max Int</b> <input type="checkbox"/>
Sport M. <b>15</b>	Relaxation, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

**Legend:** **Low** = Low - **High** = High - **Int** = Intermediate - **Max** = Maximum

# PROGRAMME LIST ELITE 2 - SPORT WOMAN

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <b>Intensity</b>
Sport W. <b>16</b>	Capillarisation	As shown in the pictures at the end of this manual	20' <b>Low</b>
Sport W. <b>17</b>	Warm-up	As shown in the pictures at the end of this manual	10' <b>Low</b>
Sport W. <b>18</b>	Recovery, post-competition/training	As shown in the pictures at the end of this manual	20' <b>Low</b>
Sport W. <b>19</b>	Maximal force, lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' <b>Int Max Int</b>
Sport W. <b>20</b>	Resistance force, lower limbs	I1-I2-I3-I4-I5-I6	5' 20' 5' <b>Int Max Int</b>
Sport W. <b>21</b>	Explosive force, lower limbs	I1-I2-I3-I4-I5-I6	5' 15' 5' <b>Int Max Int</b>
Sport W. <b>22</b>	Reactivity, lower limbs	I1-I2-I3-I4-I5-I6	5' 20' 5' <b>Int Max Int</b>
Sport W. <b>23</b>	Relaxation, lower limbs	I1-I2-I3-I4-I5-I6	20' <b>Low</b>
Sport W. <b>24</b>	Maximal force, abdominals and dorsals	T1-T2-T3-T4-T5-T6	5' 15' 5' <b>Int Max Int</b>
Sport W. <b>25</b>	Resistance force, abdominals and dorsals	T1-T2-T3-T4-T5-T6	5' 20' 5' <b>Int Max Int</b>
Sport W. <b>26</b>	Relaxation, abdominals and dorsals	T1-T2-T3-T4-T5-T6	20' <b>Low</b>
Sport W. <b>27</b>	Maximal force, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' <b>Int Max Int</b>
Sport W. <b>28</b>	Resistance force, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 20' 5' <b>Int Max Int</b>
Sport W. <b>29</b>	Explosive force, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	5' 15' 5' <b>Int Max Int</b>
Sport W. <b>30</b>	Relaxation, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' <b>Int</b>

# PROGRAMME LIST ELITE 2 - REHAB

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <span style="color: red;">Intensity</span>
Rehab <b>31</b>	DEMO	As shown in the pictures at the end of this manual	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">2' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="text-align: center;">2' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="text-align: center;">2' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>32</b>	Motor point search	As shown in the pictures at the end of this manual	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">5' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>33</b>	Analgetic Tens	As shown in the pictures at the end of this manual	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">30' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>34</b>	Endorphinic Tens	As shown in the pictures at the end of this manual	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">20' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="text-align: center;">20' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>35</b>	Contractures	As shown in the pictures at the end of this manual	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">30' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>36</b>	Hematoma	As shown in the pictures at the end of this manual	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">20' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>37</b>	Sciatic neuralgia	R7	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">30' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>38</b>	Recovery, lower limbs atrophy	I1-I2-I4	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">5' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Int</div> <div style="text-align: center;">15' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>High</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>39</b>	Backache	T3	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">20' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="text-align: center;">20' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>40</b>	Recovery, upper limbs atrophy	S1-S2-S3-S4-S5-S6-S7	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">5' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Int</div> <div style="text-align: center;">15' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Int</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>41</b>	Cervical pains	R1	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">30' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>42</b>	Periarthritis	R2	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">15' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="text-align: center;">10' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>
Rehab <b>43</b>	Epicondylitis	R3-R4	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">15' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="text-align: center;">10' <div style="width: 20px; height: 10px; background-color: black; margin: 2px 0;"></div>Low</div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> <div style="border: 1px solid black; width: 20px; height: 15px; margin-left: 10px;"></div> </div>

**Legend:** Low = Low - High = High - Int = Intermediate - Max = Maximum

# PROGRAMME LIST ELITE 2 - FITNESS MAN

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) Intensity
Fitness M. <b>44</b>	Relaxing massage <b>Bio Pulse</b>	As shown in the pictures at the end of this manual	25' Low <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>45</b>	Toning massage	As shown in the pictures at the end of this manual	20' Low <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>46</b>	Vascularisation, lower limbs	11-12-13-14-15-16	25' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>47</b>	Drainage, lower limbs	11-12-13-14-15-16	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>48</b>	Lipolysis, lower limbs	11-12-13-14-15-16	30' 20' Int Int <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>49</b>	Firming, lower limbs	11-12-13-14-15-16	20' High <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>50</b>	Toning, lower limbs	11-12-13-14-15-16	15' High <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>51</b>	Shaping and sculpting, lower limbs	11-12-13-14-15-16	20' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>52</b>	Mass development, lower limbs	11-12-13-14-15-16	15' Max <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>53</b>	Drainage, abdominals and dorsals	T1-T2-T3-T4-T5-T6	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>54</b>	Lipolysis, abdominals and dorsals	T1-T2-T3-T4-T5-T6	30' 20' Int Int <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>55</b>	Firming, abdominals and dorsals	T1-T2-T3-T4-T5-T6	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>56</b>	Toning, abdominals and dorsals	T1-T2-T3-T4-T5-T6	15' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness M. <b>57</b>	Sculpting, abdominals and dorsals	T1-T2-T3-T4-T5-T6	20' Int <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

## PROGRAMME LIST ELITE 2 - FITNESS MAN

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <b>Intensity</b>
Fitness M. <b>58</b>	Drainage, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' <div style="display: flex; align-items: center;"> <div style="width: 20px; height: 10px; background-color: black; margin-right: 5px;"></div> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px;" type="text"/> </div> <b>Int</b>
Fitness M. <b>59</b>	Lipolysis, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	30'      20' <div style="display: flex; align-items: center;"> <div style="width: 20px; height: 10px; background-color: black; margin-right: 5px;"></div> <div style="width: 20px; height: 10px; background-color: black; margin-right: 5px;"></div> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px;" type="text"/> </div> <b>Int      Int</b>
Fitness M. <b>60</b>	Firming, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' <div style="display: flex; align-items: center;"> <div style="width: 20px; height: 10px; background-color: black; margin-right: 5px;"></div> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px;" type="text"/> </div> <b>High</b>
Fitness M. <b>61</b>	Toning, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	15' <div style="display: flex; align-items: center;"> <div style="width: 20px; height: 10px; background-color: black; margin-right: 5px;"></div> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px;" type="text"/> </div> <b>High</b>
Fitness M. <b>62</b>	Sculpting, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	20' <div style="display: flex; align-items: center;"> <div style="width: 20px; height: 10px; background-color: black; margin-right: 5px;"></div> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px;" type="text"/> </div> <b>Max</b>
Fitness M. <b>63</b>	Mass development, upper limbs and pectorals	S1-S2-S3-S4-S5-S6-S7-S8	15' <div style="display: flex; align-items: center;"> <div style="width: 20px; height: 10px; background-color: black; margin-right: 5px;"></div> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px; margin-right: 5px;" type="text"/> <input style="width: 40px; height: 20px;" type="text"/> </div> <b>Max</b>

**Legend:** **Low** = Low - **High** = High - **Int** = Intermediate - **Max** = Maximum

# PROGRAMME LIST ELITE 2 - FITNESS WOMAN

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <b>Intensity</b>
Fitness W. <b>64</b>	Relaxing massage <b>Bio Pulse</b>	As shown in the pictures at the end of this manual	25' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>65</b>	Toning massage	As shown in the pictures at the end of this manual	20' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>66</b>	Vascularisation,	I1-I2-I3-I4-I5	25' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>67</b>	Drainage, hips and thighs 1	I1-I2-I3-I4-I5	20' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>68</b>	Drainage, hips and thighs 2 <b>Bio Pulse</b>	I1-I2-I3-I4-I5	20' <b>Max</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>69</b>	Lipolysis, hips and thighs 1	I1-I2-I3-I4-I5	30' 20' <b>Int Int</b> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>70</b>	Lipolysis, hips and thighs 2 <b>Bio Pulse</b>	I1-I2-I3-I4-I5	20' 20' <b>Int Int</b> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>71</b>	Firming, hips and thighs 1	I1-I2-I3-I4-I5	20' <b>High</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>72</b>	Firming, hips and thighs 2 <b>Bio Pulse</b>	I1-I2-I3-I4-I5	20' <b>High</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>73</b>	Toning, hips and thighs	I1-I2-I3-I4-I5	15' <b>High</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>74</b>	Sculpting, hips and thighs	I1-I2-I3-I4-I5	20' <b>High</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>75</b>	Drainage, gluteus 1	16	20' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>76</b>	Drainage, gluteus 2 <b>Bio Pulse</b>	16	20' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>77</b>	Lipolysis, gluteus 1	16	30' 20' <b>Int Int</b> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>78</b>	Lipolysis, gluteus 2 <b>Bio Pulse</b>	16	20' 20' <b>Int Int</b> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>79</b>	Firming, gluteus 1	16	20' <b>High</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>80</b>	Firming, gluteus 2 <b>Bio Pulse</b>	16	20' <b>High</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

# PROGRAMME LIST ELITE 2 - FITNESS WOMAN

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <span style="color: red; font-weight: bold;">Intensity</span>
Fitness W. <b>81</b>	Toning, gluteus 1	I6	15' <span style="background-color: black; color: red;">High</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>82</b>	Toning, gluteus 2	I6	15' <span style="background-color: black; color: red;">High</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>83</b>	Sculpting, gluteus 1	I6	20' <span style="background-color: black; color: red;">High</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>84</b>	Sculpting, gluteus 2	I6	20' <span style="background-color: black; color: red;">High</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>85</b>	Drainage, abdominal muscles 1	T5-T6	20' <span style="background-color: black; color: red;">Int</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>86</b>	Drainage, abdominal muscles 2 <b>Bio Pulse</b>	T5-T6	20' <span style="background-color: black; color: red;">Int</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>87</b>	Lipolysis, abdominal muscles 1	T5-T6	30' <span style="background-color: black; color: red;">Int</span> 20' <span style="background-color: black; color: red;">Int</span> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>88</b>	Lipolysis, abdominal muscles 2 <b>Bio Pulse</b>	T5-T6	20' <span style="background-color: black; color: red;">Int</span> 20' <span style="background-color: black; color: red;">Int</span> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>89</b>	Firming, abdominal muscles 1	T5-T6	20' <span style="background-color: black; color: red;">High</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>90</b>	Firming, abdominal muscles 2 <b>Bio Pulse</b>	T5-T6	20' <span style="background-color: black; color: red;">High</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>91</b>	Toning, abdominal muscles	T5-T6	15' <span style="background-color: black; color: red;">High</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>92</b>	Sculpting, abdominal muscles	T5-T6	20' <span style="background-color: black; color: red;">High</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>93</b>	Firming, breast	S9	15' <span style="background-color: black; color: red;">Low</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Fitness W. <b>94</b>	Sculpting, breast	S9	10' <span style="background-color: black; color: red;">Low</span> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

**Legend:** Low = Low - High = High - Int = Intermediate - Max = Maximum

# PROGRAMME LIST ELITE 2 - FITNESS WOMAN


Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <b>Intensity</b>
Fitness W. <b>95</b>	Drainage, upper limbs	S1-S2-S3-S4-S5-S6-S7	20' ■ Int <input type="text"/> <input type="text"/> <input type="text"/>
Fitness W. <b>96</b>	Lipolysis, upper limbs	S1-S2-S3-S4-S5-S6-S7	30' ■ Int <input type="text"/> <input type="text"/> <input type="text"/>
Fitness W. <b>97</b>	Capillarisation face	V1	15' ■ Low <input type="text"/> <input type="text"/> <input type="text"/>
Fitness W. <b>98</b>	Lifting effect, face	V1	15' ■ Low <input type="text"/> <input type="text"/> <input type="text"/>
Fitness W. <b>99</b>	Swollen legs	I1-I2-I4	15' ■ Low <input type="text"/> <input type="text"/> <input type="text"/>
Fitness W. <b>100</b>	Cramp prevention, lower limbs	I1-I2-I4	15' ■ Low <input type="text"/> <input type="text"/> <input type="text"/>

# PROGRAMME LIST GENESY 2

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <b>Intensity</b>
<b>1</b>	Demo	As shown in the pictures at the end of this manual	2'    2'    2' <input type="checkbox"/> <b>Low</b> <b>Low</b> <b>Low</b>
<b>2</b>	AnalgeticTens	As shown in the pictures at the end of this manual	30' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>3</b>	Endorphinic Tens	As shown in the pictures at the end of this manual	20'    20' <input type="checkbox"/> <input type="checkbox"/> <b>Low</b> <b>Low</b>
<b>4</b>	Tens Burst	As shown in the pictures at the end of this manual	20' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Int</b>
<b>5</b>	Analgetic Tens acute inflammation	As shown in the pictures at the end of this manual	30' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>6</b>	Analgetic Tens chronic inflammation	As shown in the pictures at the end of this manual	30' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>7</b>	Patellar chondritis (knee)	R8	30' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Int</b>
<b>8</b>	Hematoma	As shown in the pictures at the end of this manual	20' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>9</b>	Muscular contractures	As shown in the pictures at the end of this manual	20' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>10</b>	Cervical pains	R1	30' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>11</b>	Torticollis	R1	30' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>12</b>	Low backache	T3-R7	10'    10' <input type="checkbox"/> <input type="checkbox"/> <b>Low</b> <b>Low</b>
<b>13</b>	Epicondylitis	R3-R4	15'    10' <input type="checkbox"/> <input type="checkbox"/> <b>Low</b> <b>Low</b>
<b>14</b>	Scapulo-humeral periarthritis (shoulder)	R2	15'    10' <input type="checkbox"/> <input type="checkbox"/> <b>Low</b> <b>Low</b>
<b>15</b>	Carpal tunnel syndrome	Place an electrode on the front part of the sore wrist and the other electrode on the back of the hand. The polarity is irrelevant.	15' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>16</b>	Cervico-brachial neuralgia	R5	20' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>
<b>17</b>	Antistress	Place an electrode under the sole of the right foot and the other electrode under the sole of the left foot. The polarity is irrelevant.	20' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Int</b>

**Legend:** **Low** = Low - **High** = High - **Int** = Intermediate - **Max** = Maximum

# PROGRAMME LIST GENESY 2

Programme No.	Programme Name	Electrode Positioning	Phase times (min.)	
			Intensity	
18	Nausea	Place an electrode on the front part of the wrist and the other electrode on the back of the same hand. The polarity is irrelevant.	10' Low	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
19	Hip arthrosis	I3	30' Low	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
20	Menstrual pain	T6	30' Int	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
21	Muscular tone recovery, upper limbs 1	S1-S2-S3-S4	5' 15' Int Max	<input type="checkbox"/> <input type="checkbox"/>
22	Muscular tone recovery, upper limbs 2	S1-S2-S3-S4	5' 15' Int Max	<input type="checkbox"/> <input type="checkbox"/>
23	Mobilization, upper limbs	S1-S2-S3-S4	15' Int	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
24	Re-musclation upper limbs	S1-S2-S3-S4	15' Int	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
25	Muscular tone recovery, lower limbs 1	I1-I2-I3-I4-I5	5' 15' Int Max	<input type="checkbox"/> <input type="checkbox"/>
26	Muscular tone recovery, lower limbs 2	I1-I2-I3-I4-I5	5' 15' Int Max	<input type="checkbox"/> <input type="checkbox"/>
27	Mobilization, lower limbs	I1-I2-I3-I4-I5	15' Int	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
28	Re-musclation lower limbs	I1-I2-I3-I4-I5	15' Int	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
29	Muscular tone recovery, shoulders	S5-S6-S7	5' 15' Int Max	<input type="checkbox"/> <input type="checkbox"/>
30	Mobilization, shoulders	S5-S6-S7	15' Int	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
31	Re-musclation shoulders	S5-S6-S7	15' Int	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
32	Muscular tone recovery, ankle		5' 15' Int Max	<input type="checkbox"/> <input type="checkbox"/>
33	Mobilization, ankle		15' High	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
34	Re-musclation ankle		15' High	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

# PROGRAMME LIST GENESY 2

Programme No.	Programme Name	Electrode Positioning	Phase times (min.) <b>Intensity</b>
<b>35</b>	Hemiplegia	S1-S2-S3-S4-S6-S7-I1-I2-I4	20' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>36</b>	Spasticity	S1-S2-S3-S4-S6-S7-I1-I2-I4	20' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>37</b>	Muscle strengthening, foot (flatfoot)	Place the electrodes under the sole of the foot	15' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>38</b>	Incontinence, Mixed	Use with probes in compliance with the law	15' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>39</b>	Incontinence, Vesical	Use with probes in compliance with the law	15' <b>Low</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>40</b>	Incontinence, from stress	Use with probes in compliance with the law	15' <b>Int</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>41</b>	IONOPHORESIS	Place sponge electrodes respecting the correct polarity (as indicated on drug packaging)	20' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>42</b>	Maximal force,	As shown in the pictures at the end of this manual	5'    15'    5' <input type="checkbox"/> <b>Int</b> <b>Max</b> <b>Int</b>
<b>43</b>	Resistance force,	As shown in the pictures at the end of this manual	5'    15'    5' <input type="checkbox"/> <b>Int</b> <b>Max</b> <b>Int</b>
<b>44</b>	Mass development	As shown in the pictures at the end of this manual	5'    15'    5' <input type="checkbox"/> <b>Int</b> <b>Max</b> <b>Int</b>
<b>45</b>	Vascularisation,	As shown in the pictures at the end of this manual	20' <b>High</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>46</b>	Drainage,	As shown in the pictures at the end of this manual	20' <b>High</b> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<b>47</b>	Lipolysis,	As shown in the pictures at the end of this manual	20'    20' <input type="checkbox"/> <input type="checkbox"/> <b>High</b> <b>High</b>
<b>48</b>	Firming,	As shown in the pictures at the end of this manual	5'    15' <input type="checkbox"/> <input type="checkbox"/> <b>High</b> <b>High</b>
<b>49</b>	Toning,	As shown in the pictures at the end of this manual	5'    20' <input type="checkbox"/> <input type="checkbox"/> <b>High</b> <b>High</b>
<b>50</b>	Cramp prevention	As shown in the pictures at the end of this manual	15' <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <b>Low</b>

**Legend:** **Low** = Low - **High** = High - **Int** = Intermediate - **Max** = Maximum

## PROGRAMME LIST REHAB 2

Programme No.	Programme Name	Electrode Positioning	Phase times (min.)			
			<b>Intensity</b>			
<b>1</b>	Endorphinic Tens	As shown in the pictures at the end of this manual	20' Low	20' Low	<input type="checkbox"/>	<input type="checkbox"/>
<b>2</b>	Analgetic Tens	As shown in the pictures at the end of this manual	30' Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>3</b>	Cervical pain	R1	30' Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>4</b>	Cervico-brachial pain	R5	15' Low	10' Low	<input type="checkbox"/>	<input type="checkbox"/>
<b>5</b>	Backache	T3	10' Low	10' Low	<input type="checkbox"/>	<input type="checkbox"/>
<b>6</b>	Sciatic neuralgia	R7	30' Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>7</b>	Post-accident muscle strengthening, Ant. Cruciate Ligament	I1-I2-I3-I4-I5	5' High	15' Int	5' High	<input type="checkbox"/>
<b>8</b>	Atrophy, quadriceps	I4	5' High	15' Int	5' High	<input type="checkbox"/>
<b>9</b>	Atrophy, gemellus	I1	5' High	15' Int	5' High	<input type="checkbox"/>
<b>10</b>	Atrophy, shoulder	S5-S6-S7	5' High	15' Int	5' High	<input type="checkbox"/>

## PROGRAMME LIST REHAB 2

Programme No.	Programme Name	Electrode Positioning	Phase times (min.)			
			<b>Intensity</b>			
<b>11</b>	Atrophy, biceps	As shown in the pictures at the end of this manual	5' High	15' Int	5' High	<input type="checkbox"/>
<b>12</b>	Muscle strengthening	As shown in the pictures at the end of this manual	5' High	15' Int	5' High	<input type="checkbox"/>
<b>13</b>	Relaxation,	As shown in the pictures at the end of this manual	20' Int	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>14</b>	Cramp prevention, lower limbs	As shown in the pictures at the end of this manual	15' Int	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>15</b>	Incontinence, from stress	As shown in the pictures at the end of this manual	15' Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>16</b>	Stress incontinence	As shown in the pictures at the end of this manual	15' Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>17</b>	Incontinence, Mixed	I1-I2-I3-I4-I5	15' Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>18</b>	Limb swelling	As shown in the pictures at the end of this manual	15' Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>19</b>	Lipolysis,	As shown in the pictures at the end of this manual	20' High	20' High	<input type="checkbox"/>	<input type="checkbox"/>
<b>20</b>	Firming,	I1-I2-I3-I4-I5	5' High	20' High	<input type="checkbox"/>	<input type="checkbox"/>
<b>21</b>	IONOPHORESIS	Place sponge electrodes respecting the correct polarity (as indicated on drug packaging)	20'	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### UROLOGY PROGRAMMES INFORMATION (Probes are available upon payment)

#### Type

Urology applications involve the use of endovaginal or endorectal probes for the specific intended use and which must be EC MDD certified according to the 93/42/EEC directive; probes must be bipolar and equipped with a 2 mm female socket for cables with 2 mm male connectors.

#### Warnings

Given the medical use of these applications, they should only be performed upon consent of medical personnel.

#### Use

For a correct use of probes, follow the instructions provided by the manufacturer or medical personnel supervising the patient.

#### Maintenance

**For cleaning, sterilizing, and disinfecting the probes refer to the manufacturer's instructions.**

#### Recommendations

In case of deterioration of the probe, it should be replaced and no longer used.

# FREQUENTLY ASKED QUESTIONS

*Why does the intensity return immediately to 0 when the value is increased to over 10 mA?*

As already described in this manual, the device is provided with a current supply control system. If the intensity value is increased to more than 10 mA and the circuit is open, the electrostimulator brings the intensity value back to 0 mA. Thus it is necessary, before beginning any treatment, to make sure that the cables are connected to the device, that the electrodes are applied to the area to be treated and that they are in good condition so as not to prejudice their conduction capacity

*Which electrodes should be used?*

It is recommended to use self-adhesive electrodes that are more handy and improve stimulation quality. If used properly (e.g., on clean skin), the electrodes can be used for 25-30 applications. They should, however, always be replaced if they are not perfectly adherent to the skin.

*Where should the electrodes be placed?*

This manual includes pictures related to electrode positioning for all body parts (polarity is irrelevant). It is thus sufficient to follow these indications.

The correct positioning may be verified using the special motor point pen or the following method: apply the electrodes as indicated in the picture and start the stimulation, then use one hand to push the electrode in various directions making the skin slide over the muscle. You will notice an increase or decrease in the stimulation depending on the position of the electrode. Once the point of greatest stimulation has been identified, lower the channel intensity to zero (0.0 mA), replace the electrode and gradually increase the intensity.

*Use of split cables and Y cables*

Y-cables double the number of electrodes per channel. This allows to work, for example, on the quadriceps vastus medialis and lateralis with the same channel, thus using both split channels to perform the treatment simultaneously on two limbs involving four muscles. Use for medical applications is not recommended.

*Does the use of Y-cables entail a decrease in power?*

The current intensity supplied for each channel does not vary. However, when using Y cables on a channel, the current is distributed over a larger muscle mass and the contraction reduced. To obtain the same contraction, increase the intensity.

*Can the electrostimulator be harmful?*

The device is unlikely to cause muscle damage. An important principle to follow during stimulation is to gradually increase the intensity while paying attention to the behaviour of the muscle, avoiding complete limb extension. In case of doubt, consult a specialist.

*Is it possible to use the electrostimulator during the menstrual cycle?*

Interferences such as early or late periods or heavier or lighter flows are extremely subjective and variable. It is however recommended to avoid treatments in the abdominal region during and immediately before-after the cycle.

*Is it possible to use the electrostimulator while breast-feeding?*

As of today no side-effects have been reported concerning the use of electrostimulation while breast-feeding. It is recommended, however, not to treat the thoracic region.

*Does the presence of dermatological diseases such as psoriasis, urticaria, etc. prevent the use of the electrostimulator?*

Yes, it is definitely advisable not to treat epidermic regions with severe dermatological pathologies.

*How long does it take for the first results to become apparent?*

Cosmetic results are clearly subjective. It is however possible to state that good toning results may be achieved after just 15 days with 3 or 4 regular sessions per week; whereas electrolipolysis and electrodrainage treatments require 40 days. Better and quicker results can be achieved if the treatments are combined with good physical activity and a healthy lifestyle.

*How many weekly electrostimulation sessions can I perform?*

Regarding physical preparation sessions, we recommend that you refer to the weekly programmes included in the Personal Trainer. The weekly frequency of fitness and beauty applications depends on the type of treatment. Toning treatments should be performed 3-4 times a week on alternate days. Lipolysis and drainage treatments can instead be performed on a daily basis.

# GUARANTEE

The electrostimulator is guaranteed to the first user for 24 months from the purchase date against material or manufacturing defects, provided that it is used properly and maintained under normal operating conditions.

In order to take advantage of the guarantee service, the user must comply with the following guarantee clauses:

1. The products to be repaired must be sent by and at the expenses of the customer in their original packages.
2. The product's guarantee is subject to the production of a fiscal document (fiscal receipt, receipted bill or sale invoice), attesting to the product's purchase date.
3. The repair work shall have no effect on the original guarantee expiry date and shall neither renew nor extend it.
4. If, prior to carrying out the repair work, no defects are found, the costs related to inspection times shall in any case be charged.
5. The guarantee lapses if the breakdown has been caused by: bumps, falls, incorrect or improper use of the product, use of an external non-original power supply, accidental events, modification, replacement/detachment of guarantee seals and/or product tampering. Said guarantee does not cover damage caused during transport when unfit packages are used (see point 1).
6. The guarantee does not cover the inability to use the product, other incidental or consequent costs or other expenses incurred by the purchaser.

Note: Before returning the device for repairs, read carefully the user instructions contained in the manual. In case of any damage caused to the device, refer to your dealer. The manufacturer reserves the right to make all necessary modifications at any time in order to improve the appearance and quality of the product.

Manufactured by:

DOMINO S.R.L  
Tel. +390438793052  
Fax +390438796463  
E-Mail: info@domino.tv.it

Distributed by:

Globus Italia  
Tel. +39-043879333  
Fax +39-0438793363  
E-Mail: info@globusitalia.com

# Purchasable accessories

Motor Point Pen

Gel

Ionophoresis electrode

Battery charger

Y cables

Split cables

Urology probe

Urology probe

*The size and features described in this manual may differ and are not binding. The manufacturer reserves the right to make modifications without notice.*

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