SLENDERTONE®

United Kingdom 0845 070 77 77
Republic of Ireland 1890 92 33 88
France 0810 347 450
España: 900994467
International: +353 1 844 1016

info@slendertone.com

Designed by & Manufactured for:

Bio-Medical Research Ltd., Parkmore Business Park West, Galway, Ireland.
Copyright: © 2005, Bio-Medical Research Ltd. All Rights Reserved
Part No: 2400-5401 Rev.: 10 Date of Issue: 11/05

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www.slendertone.com

If this is your first ACTIVE session, you must read the rest of this manual before you begin.

WHAT IS SLENDERTONE ACTIVE?

SLENDERTONE ACTIVE is a breakthrough in personalised health and fitness technology, which tracks your heart rate to ensure you exercise at the best pace to burn fat and improve fitness. As you exercise, ACTIVE delivers a specially designed muscle toning treatment, using an abdominal belt, which strengthens, firms and flattens your whole abdominal area while you exercise.

THE SLENDERTONE ACTIVE THEORY

We all need exercise in order to stay fit and healthy. However, how do you know whether or not you are exercising at the right level? One tried and tested method is to measure your heart rate and this is exactly what SLENDERTONE ACTIVE does for you. It monitors your heart rate as you exercise. You can then check your display to see what your heart rate is and you'll know whether you are exercising below, above or at your optimum level.

Using this technique, ACTIVE ensures that you always exercise at peak efficiency. The reason for this is simple. When you exercise with your heart rate at its optimum range (your "Training Zone"), your fitness improves.

QUICK GUIDE TO USING SLENDERTONE ACTIVE

As an example; if you wish to use the Walk and Tone plan (P1), just follow the steps below.

(a) SET-UP

1. Setup and switch on your ACTIVE unit (see page 6).
2. Increase the intensity from the belt to a level you find comfortable.

(b) GO

3. Go! Walk as the unit tones your abs.
4. Complete your session, checking your heart rate on the display from time to time. Make sure you are exercising at the correct level – see the heart rate table. e.g. if you are aged 45 keep your heart rate between 106 and 141.

(c) STOP

5. Once your session is over, complete the training diary with your ACTIVE points score for that session. The ACTIVE points table tells you how many points to fill in.

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6. Reach your recommended weekly points target and feel the results! In the example below the target is 20 ACTIVE points. We recommend this target for beginners.

SLENDERTONE ACTIVE Diary / Journal / Trainingsplan / Diario / Dagboek											
Day / Jour / Tag / Día / Dag / Giorno											
1	2	3	4	5	6	7					
4/10/06	5/10/06	6/10/06	7/10/06	8/10/06	9/10/06	10/10/06					
4	6	6	6	4							
											20 ★

Age	Walk & Tone
20 - 29	114 - 152 ♥
30 - 39	110 - 147 ♥
40 - 49	106 - 141 ♥
50 - 59	102 - 136 ♥
60 - 69	98 - 130 ♥

HEART RATE TABLE

Time	Walk & Tone
15:00	2 ☆
30:00	4 ☆
45:00	6 ☆
60:00	8 ☆

ACTIVE POINTS TABLE
(FOR WALK & TONE)

Replacement SLENDERTONE ACTIVE diaries are available on www.slendertone.com.

DETERMINE YOUR READINESS FOR PHYSICAL ACTIVITY

As with any exercise, before starting a SLENDERTONE ACTIVE training plan you must first ensure you are physically ready for exercise. If you are in any doubt, please check the contra-indications on page 11 and if necessary, consult your doctor before beginning.

Please answer the questions below.

1. Has your doctor ever warned you of a heart condition or said that you should only do exercises recommended by a doctor? ☐ Yes / No ☐
2. Do you feel chest pains when you exercise? ☐ Yes/ No ☐
3. In the past month, have you experienced any chest pain while you were not participating in any physical activity? ☐ Yes/ No ☐
4. Do you lose your balance because of dizziness or do you ever lose consciousness? ☐ Yes/ No ☐
5. Do you have any bone or joint problems which could be made worse by a change in your level of physical activity? ☐ Yes/ No ☐
6. Are you currently taking medication for blood pressure or for an underlying heart condition? ☐ Yes/ No ☐
7. Do you know of any other reason why you should not participate in physical activity? ☐ Yes/ No ☐

If you answer yes to one or more of the above questions, you should consult your doctor before starting a SLENDERTONE ACTIVE training plan. SLENDERTONE recommends that you consult your doctor before undertaking the physical exercises mentioned in this manual. As with any new exercise, remember to always begin slowly and gradually build-up the level of exercise.

YOUR SLENDERTONE ACTIVE PACK

1. ACTIVE Belt (fits waist 28" – 42"):

With the three SLENDERTONE pads (1a) positioned on the belt (page 6), fasten the belt securely but comfortably around your waist.

2. Base Unit:

Slide the cover off the base unit and insert the batteries. Then close the cover and slide the base unit into the holster on the belt.

3. Heart Rate Pad:

The heart rate pad is placed on the chest strap to monitor your heart rate during each SLENDERTONE ACTIVE session.

4. Heart Rate Lead:


This lead is used to connect the belt to the chest strap.

5. Chest Strap (fits chest size 24" – 47"):

Fasten this strap around your chest in order to hold the pad firmly in place.

6. Control Unit:

This is the brain of your SLENDERTONE ACTIVE unit. It controls all of the functions described overleaf.

SLENDERTONE ACTIVE CONTROL UNIT

ON/ Off Button (⏻):

Press for 2 seconds to switch your unit on or off. Press briefly during a session to pause the session and press briefly again to restart it.

Personalisation Button (⚙):

Press this button to activate the personalisation screen. Use either intensity button to scroll to the **code** which best corresponds to you and then press the personalisation button again to enter.

Example:

If John is 44 years old (read line across from 40-49) and weighs 65 Kg (read down from 60-69), you can determine, using this chart, that John's **code** is 9. Now determine your **code** from the chart.


Code Table:

		Weight (Kg)				
		50 - 59	60 - 69	70 - 79	80 - 89	90 +
Age	20 - 29	15	15	14	13	13
	30 - 39	12	12	11	10	10
	40 - 49	9	9	8	8	7
	50 - 59	6	5	5	4	4
	60 - 65	3	3	2	2	1

Personalising gives extra heart rate feedback and improves the accuracy of your calorie counter.

NOTE:

Calorie feedback is an estimate and is based on keeping your heart rate within your Training Zone.

Information Button (i):

Press this button to check your average heart rate, calories burned, ACTIVE points, intensity levels and total number of sessions (see page 7).

Training Plan Button (P):

Select one of the three Training Plans using this button (see page 5).

Intensity Buttons (▲/▼):

Press the up/down arrows to increase/decrease the muscle training intensity.

SLENDERTONE ACTIVE DISPLAY

The display on your SLENDERTONE ACTIVE control unit supplies you with all the information you need to exercise at your ideal pace. This will let you keep track of your progress.

- Graphical heart rate indicator.
- Heart rate Training Zone.
- 123 ♥ Your current heart rate.
- Avg 123 ♥ Average heart rate in current session.
- 10:08 Time remaining/ elapsed in current session.
- 150 Intensity on your left side.
- 46 Intensity on your right side.
- T Highest intensity for current/ last session.

- kcal 123 Number of calories burned in the current session.
- ★45 ACTIVE points scored in current session.
- || Programme is paused.
- Lock Keylock function activated.
- Mute Mute function activated.
- Battery Battery power remaining.
- # Number of sessions to date.
- △ Poor pad contact.

SLENDERTONE ACTIVE TRAINING PLANS

SLENDERTONE ACTIVE offers you a choice of three training plans. Each plan will benefit you in a different way. For example, if you want to improve your health and fitness, you can use either Plan 1 or Plan 2, depending on your current level of fitness. On the other hand, if you want to just tone your abdominal muscles, you can use Plan 3.

Each plan has a weekly points target. When you use your SLENDERTONE ACTIVE you earn ACTIVE points. The table below shows you the weekly targets for each of the 3 training plans.

NOTE:

SLENDERTONE ACTIVE may also be used with a treadmill or a stationary exercise bike.

Plan	Description	What to do	Weekly target
P1	Walk & Tone	Go for a brisk walk while using the SLENDERTONE ACTIVE belt 3-5 times per week.	20 ★
P2	Jog & Tone	Go jogging while using the SLENDERTONE ACTIVE belt 2-4 times per week.	18 ★
P3	Tone	Use the SLENDERTONE ACTIVE belt at work or in the comfort of your home 5 times per week for 30 minutes.	5 ★

SLENDERTONE ACTIVE POINTS

The **ACTIVE Points Table** below shows you how many points you get for each of your SLENDERTONE ACTIVE Training Plans. You will also receive more points by exercising for longer within the individual plans. As you can see from the table below, 30 minutes of Walk & Tone earns you 4 points, whereas 45 minutes of Walk & Tone earns you 6 points. Feel free to go over your weekly target if you wish to get even more from your SLENDERTONE ACTIVE unit.

When exercising, you must keep your heart rate at the correct level. The **Heart Rate Table** below shows you what your heart rate should be when using the Walk & Tone and Jog & Tone plans, depending on your age. This is called your Training Zone. This is necessary to train at the best pace to burn fat. For best results, use the highest muscle training intensity you can.

ACTIVE POINTS TABLE (★)


MINUTES EXERCISED	WALK & TONE	JOG & TONE	TONE
15:00	2 ★	3 ★	1/2 ★
30:00	4 ★	6 ★	1 ★
45:00	6 ★	9 ★	- ★
60:00	8 ★	12 ★	- ★
1 ★ ≈ 50 kcal.			

HEART RATE TABLE (♥)

Age	Walk & Tone	Jog & Tone
20 - 29	114 - 152 ♥	143 - 168 ♥
30 - 39	110 - 147 ♥	138 - 161 ♥
40 - 49	106 - 141 ♥	132 - 155 ♥
50 - 59	102 - 136 ♥	127 - 149 ♥
60 - 69	98 - 130 ♥	122 - 143 ♥

USING YOUR ACTIVE UNIT

1. Remove the cover from the **patterned** side of the large adhesive pad (Fig. a). Place the pad between the parallel lines in the centre of the belt. Ensure the pad completely covers the metal stud.
2. Remove the covers from the **patterned** side of the 2 smaller pads. Position these either side of the large pad within the markings best suited to your figure (Fig. b). Make sure the studs are fully covered by the pads.

3. Remove the covers from the **black** side of the 3 pads. Position the belt with the large pad over your belly-button and the small pads between your hip-bones and ribs and fasten the belt securely but comfortably around your waist.

NOTES:

Do not throw any of the covers away, as you will need them later. Ensure you position the belt the right way up.

4. Open the battery cover on the base unit and insert 3 AAA batteries, as indicated in the base of the compartment (Fig. c). Push the base unit into the holder until the guidelines are correctly aligned. This may require some force (Fig. d).

5. Attach heart rate lead to the chest strap (fig. e) and close the outer cover over it (Fig. f).

6. Remove the cover from the patterned side of the heart rate pad and place it on the inside of the chest strap, directly over the metal contact (Fig. g). Remove the cover from the black side of the heart rate pad (do not throw any of the covers away) and fasten the chest strap around your chest with the pad just under your heart (Fig. h). Adjust the chest strap to fit comfortably around your chest and secure it using the buckle. You may need to try this a couple of times in order to get the correct position for a strong heart rate signal.

NOTE:

If you have any problems, adjust the position of the chest strap to get better contact with your chest.

7. Press the on/ off button for 2 seconds to switch your unit on.

8. Select your required training plan using the training plan button (P) and then increase the intensity until you feel your abs being exercised.

NOTE:

You do not have to reach the same intensity on both sides of your body.

PAUSING THE SESSION:

Press the on/ off button briefly if you wish to pause the session. Press it briefly a second time to restart.

Fig. a


Fig. b


Fig. c

Fig. d


Fig. e


Fig. f


Fig. g

Fig. h

9. At the end of a session, the control unit's display will alternate between Fig. a and Fig. b. You can use this information to record your average heart rate for that session, the number of calories you have burned and the ACTIVE points earned. You can also check your heart rate at the end of the session and whether or not you were exercising within your training zone.

THE INFORMATION BUTTON (i)

Your SLENDERTONE ACTIVE unit stores information about current and previous sessions which can be accessed by using the information button. You can access this information while a programme is running.

Press the information button once to see (Fig. c):

1. Your average heart rate for the current session.
2. The number of calories you've burnt during the current session.
3. Your ACTIVE points score for the current session.

Press the information button twice to see (Fig. d):

4. The total number of sessions you have completed (returns to zero after 199).
5. The highest intensity you have reached in the current session.

Press the information button a third time to return to the default display.

OTHER FUNCTIONS

Intensity Keylock Function

If you find a comfortable intensity level, press the training plan button and then press the information button within 3 seconds to activate the keylock function. Do the same again to deactivate this function.

Mute Function

If you want to switch off your unit's sound effects, press the personalisation button for 2 seconds. The  symbol will appear on the display. Do the same again to deactivate the mute function.

ADDITIONAL INFORMATION

THE BELT CLIP

Your SLENDERTONE ACTIVE unit is supplied with a magnetic belt clip. Simply press the rounded side of the clip onto the back of control unit (Fig. e) and then attach the belt clip to your SLENDERTONE ACTIVE belt or another item of clothing.

ADHESIVE PADS

As well as the heart rate pad for the chest strap, you will also find 2 other heart rate pads in the pack. These pads can be used on individuals with no body hair in the upper chest area. To use one of these instead of the chest strap and pad;

1. Attach the heart rate pad to the lead (Fig. f).
2. Remove the cover from the adhesive side of the heart rate pad (Fig. g) and position the pad on the upper-right side of your chest bone. You may need to try this a couple of times in order to get the correct position for a strong heart rate signal.

Fig. a


Fig. b


Fig. c


Fig. d


Fig. e

Fig. f



Fig. g

CARING FOR YOUR BELT

You should not let your unit get wet, but you may wipe it clean with a lightly dampened cloth. The belt can be washed, but you must first remove the unit and pads. Always follow the instructions on the label when washing the belt.

Never machine wash your belt. Always hand wash it in lukewarm water. This protects the internal wiring from damage in the washing machine. You should take care even when hand washing and never wring the belt to remove water.

Do not use any bleach when washing your belt.

Do not dry clean your belt.

Do not tumble dry your belt. Dry the belt on a flat surface, but not on anything hot (e.g. a radiator) as the belt contains plastic parts. Ensure the belt is completely dry before use.

The belt must not be ironed.

CARING FOR YOUR CHEST STRAP

Wipe with a lightly dampened cloth regularly to remove any dirt/ residue.

QUESTIONS AND ANSWERS

How does ACTIVE train the abdominal muscles?

- (i) ACTIVE sends gentle signals between the pads on the belt.
- (ii) These signals travel into the muscle and switch on the nerves in that area.
- (iii) The nerves, in turn, switch on the muscles and the muscles contract.

NOTE: We do not recommend more than 60 minutes of abdominal muscle training per day.

Will the product cause muscle soreness?

As with all exercise some muscle soreness can be expected. Use a low toning intensity on your Active if you are experiencing muscle soreness.

Can I use ACTIVE for post natal exercise?

Yes, but you must wait a minimum of six weeks after childbirth and consult your doctor first. If you have had a Caesarean in the past three months, consult your doctor before using ACTIVE.

The stimulation is uncomfortable. How can I improve this?

Ensure the pads are correctly positioned and pressed firmly against your skin. Tighten the belt around your waist. Remember to switch off your unit before repositioning the belt.

How do I know when to replace the pads?

If the signal is weakening, even though the batteries are still OK, it usually indicates that the pads are worn and need replacing. Check this by inserting new batteries and/or adjusting the belt to ensure correct positioning. **Replacement pads can be purchased from your nearest SLENDERTONE stockist, from www.slendertone.com or by calling your local SLENDERTONE Careline.**

What if the heart rate reading is erratic or totally absent?

Move the chest strap upwards trying positions that are more to the side or center and place on a area where there is less body hair. Check to see that the heart rate pad is securely fastened and correctly positioned. If the heart rate pad has slipped, pause your unit and adjust it. Tighten the belt around your waist. Use a cotton shirt instead of a synthetic shirt. Use a tighter shirt to avoid fluttering of the material. Make sure you are not standing too close to any electrical equipment which may interfere with the signal. Check the batteries and replace them if necessary.

QUESTIONS AND ANSWERS

What if the heart rate reading becomes extremely high (e.g. above 180)?

Your ACTIVE unit is a very sensitive piece of equipment and may be temporarily affected by strong electro-magnetic signals. This may cause erratic and often elevated heart rate displays. Stop exercising and wait 30 seconds to see if the heart rate returns to a stable level. It is also recommended that you continuously monitor your heart rate on the display. If you are exercising too far over your recommended maximum heart rate, you should stop exercising and rest until your heart rate drops to a more normal level.

Your Maximum Heart rate = 220 minus your age!

The pads won't stick to the belt even though they are fairly new. Why is this?

Ensure the *patterned* side of each pad is placed on the belt. The *black* side of the pads goes onto your skin. Press the edges of each pad very firmly onto the belt before and after each session.

My skin is red after a session. Why?

Some skin redness after a session is natural. It is partly due to an increase in blood-flow and should fade after the treatment. It can also occur due to the pressure of the stretched belt. You should not be concerned about this, but should only wear the belt during a session. If it is excessive, you may have the intensity too high. If the problem persists, you should stop using the unit.

SLENDERTONE RECOMMENDATIONS

- We recommend that you pause the stimulator while walking up steps or over uneven ground.
- For first time users, muscle stimulation can be an unusual sensation. We recommend that you familiarise yourself with muscle stimulation at rest before progressing to muscle stimulation while walking.
- For hygiene reasons this device is for single person use only. Do not share your belt with anyone else.
- Take care while reading the display or adjusting the product while moving. These actions can result in you being less aware of your surrounding environment.
- Do not connect the units to any other device, for example portable music systems.
- Do not over-exert yourself during exercise / muscle stimulation. Any workout should be at a level that is comfortable for you.
- You should leave at least six hours between sessions.

I can feel a tingling in my legs during a session. What should I do?

Move the two small pads upward and inward on your waist (i.e. towards the centre of your body).

My waist muscles are exercising but not my stomach muscles. Why?

Move the large pad down slightly. If this doesn't help, move the small pads to a smaller figure setting and/ or increase the intensity.

Can I use ACTIVE to treat muscles weakened from lack of use due to injury?

Your ACTIVE belt may additionally be used for the alleviation of or compensation for injury. Consultation with your doctor or physio is required to establish a rehabilitation programme with the ACTIVE unit, which would safely provide improvement to the strength and tone of the abdominal muscles.

What should I do if I see EooX on the display?

Eoo2 or Eoo6 appears on the display:

- The batteries are dead.
- Replace the batteries.
- If the problem persists with new batteries return your unit to SLENDERTONE

Eoo1, Eoo3, Eoo4 or Eoo5 appears on the display:

- Internal fault.
- Return your unit to SLENDERTONE for repair.

TROUBLE-SHOOTING GUIDE

PROBLEM	POSSIBLE CAUSE	SOLUTION
The display doesn't come on & there is no signal from the unit	The batteries are inserted incorrectly The batteries are dead Low-quality batteries used	Reposition the batteries as described on page 3 Replace the batteries Use only high-quality batteries
The display is on but there's no signal	The unit has been paused No intensity set	Press the on/off button briefly Increase the intensity
The battery symbol is flashing	The batteries are low	Replace the batteries
Unpleasant feeling beneath the pads	Pads not covering the studs The pads are worn Too many consecutive sessions Pads not in contact with skin	Reposition the pads Contact your local stockist for new SLENDERTONE pads You should leave at least six hours between sessions See page 6 for correct positioning of the pads & belt
The △ symbol has appeared on the display *	Unit is not fully inserted in holster Pads not covering the studs The pads are worn The pads are on the wrong way around The covers are still on the pads Belt is loose or open on body	Push the unit fully into the holster Reposition the pads Contact your local stockist for new SLENDERTONE pads Ensure the black surface of each pads is placed on your skin Ensure the covers are removed from both sides of the pads Refasten the belt properly
The contractions are very weak even when the intensity is high	The batteries are low The pads are worn Poor pad positioning Pads not covering the studs Low quality batteries used	Replace the batteries Contact your local stockist for new SLENDERTONE pads Reposition the pads Reposition the pads Use only high-quality batteries
The unit is on but there is no heart rate reading	Chest strap too loose Heart rate lead not connected	Tighten chest strap Ensure lead connected to chest strap and belt stud

* The △ symbol always appears in conjunction with the left/right arrows. If the left arrow appears, there is a problem with the left pad, the right arrow indicates a problem with the right pad and if both arrows appear, there is a problem with either the middle pad, all three pads or the unit is not properly inserted into the holster.

CONTRA-INDICATIONS

SLENDERTONE ACTIVE is suitable for use by all healthy adults. However, as with other forms of exercise, some care is needed when using it, so always follow the points below and read the support pack before use. Some of the points below are gender specific.

Please do not use if:

- You have an electronic implant (e.g. cardiac pacemaker or defibrillator) or suffer from any other heart problem.
- You are pregnant.
- You suffer from cancer, epilepsy or are under medical supervision for cognitive dysfunction.
- The unit is in close proximity (e.g. 1m) to shortwave or microwave therapy equipment.
- You are connected to high-frequency surgical equipment.
- Wearing the device necessitates placement over areas where drugs/medicines are administered by injection (short term or long term) e.g. hormone treatment.

Please wait before using your SLENDERTONE ACTIVE until:

- At least six weeks after the birth of your baby (you must consult your doctor first).
- One month after an IUD contraceptive device (e.g. coil) has been fitted.
- At least three months after having a Caesarean (you must consult your doctor first).
- The heavy days of your period have finished.

Please get your doctor's or physio's permission before using this product if:

- If you're feeling light headed/faint please slow down or stop and pause your stimulation until you feel ready to begin again. Consult your doctor if this happens.
- You wish to use a SLENDERTONE belt and have a bad back. You should also ensure the intensity is kept low.
- You have any serious illness or injury not mentioned in this guide.
- You have recently had an operation.
- You take insulin for diabetes.
- You want to use it on a young child.
- You suffer from muscle or joint problems.
- Using the device as part of a rehabilitation programme.

When applying the pads, belt and chest strap always remember to:

- Place the pads and belt ONLY as indicated in this manual.
- Avoid placing the pads or belt on the front or sides of the neck, across or through the heart (i.e. one pad on the front of the chest and one on the back), in the genital region or on the head. (Other SLENDERTONE units are available for other areas of the body - details are available from your stockist)
- Avoid any recent scars, broken or inflamed skin, areas of infection or susceptibility to acne, thrombosis or other vascular problems (e.g. varicose veins), or any parts of the body where feeling is limited.

- Avoid areas of injury or restricted movement (e.g. fractures or sprains).
- Avoid placing the pads directly over metal implants. They may be placed on the nearest muscle.

Possible adverse reactions:

- A small number of isolated skin reactions have been reported by people using muscle stimulation devices, including allergies, a prolonged reddening of the skin and acne.

To reposition pads during a session:

- Always pause the programme currently running, unfasten the belt and then refasten it behind your back once the pads' position has been adjusted.

After strenuous exercise or exertion:

- Always use a lower intensity to avoid muscle fatigue.

Contact SLENDERTONE or an authorised distributor if:

- Your unit is not working correctly. Do not use in the meantime.
- You experience irritation, skin reaction, hypersensitivity or other adverse reaction. You should, however, note that some reddening of the skin can appear under the belt during and for a short time after a session.

NOTE:

- An effective treatment should not cause undue discomfort.

Important:

- The heart rate is an estimate and does not provide diagnostic information. Consult your doctor if in doubt.
- Rain may interfere with this device while using it outdoors. Ensure that this device is kept dry at all times.
- Keep your unit out of the reach of children.
- The studs, leads and pads must not be connected to other objects.
- Do not use your unit at the same time as any other device which transfers an electrical current into the body (e.g. another muscle stimulator).
- Do not touch the pads or metal studs while the unit is switched on.
- Avoid using while driving, crossing in traffic, operating machinery or cycling
- Your SLENDERTONE ACTIVE unit should not be used with a garment from any other SLENDERTONE product. Doing so may result in an unbalanced stimulation and an inefficient toning session.
- SLENDERTONE will not accept responsibility if the guidelines and instructions supplied with this unit are not followed.

NB: If you are in any doubt about using a SLENDERTONE belt for any reason, please consult your doctor before use.

TECHNICAL SPECIFICATIONS

Caring for your unit
Your unit should not be left in excessive sunlight. Do not allow the interior of your unit to become wet. Do not use detergents, alcohol, spray aerosols or strong solvents on your unit.
Access to the unit's interior is not required for maintenance purposes.
If your unit is damaged, you should not use it but should return it to SLENDERTONE or your local distributor for replacement or repair. Repairs, service and modifications may not be carried out by anyone other than qualified service personnel authorised by SLENDERTONE.

NOTE: Using leak-proof batteries is advised. You should never leave the batteries inside your unit if you do not intend to use it for a long period of time. If you do, the batteries may leak and damage your unit. You should be aware that some batteries sold as 'leak-proof' can still release corrosive substances, which may damage your unit. Under no circumstance should anything other than the correct type of batteries 1.5V (LR03 / AAA) be used with your unit. You must only insert the batteries in the manner specified in this manual.

Disposing of the pads and batteries
Used pads and batteries must never be disposed of in a fire but in accordance with your country's national laws governing the disposal of such items.

⚠ The battery compartment must be closed when the unit is on.
When the ⚠ symbol appears the stimulus is significantly reduced.

The signal gradually increases to a peak intensity level at the start of a contraction phase (i.e. ramp up) and gradually decreases to nothing at start of a relaxation phase (i.e. ramp down).

Accessories
When ordering new pads, you should only use those carrying the SLENDERTONE brand. Others may not be compatible with your unit and could degrade the minimum safety levels.

- SLENDERTONE Adhesive Pads:**
- 1 large adhesive pad Type 706 or 709
 - 2 small adhesive pads Type 707 or 710
 - 3 heart rate pads Type 711 x 1; 712 x 2

Batteries (3 x 1.5V AAA)

SLENDERTONE ACTIVE Belt

Chest Strap

SLENDERTONE Belt Extension

Intended use: Muscle stimulator/ Heart Rate Monitor

Waveform: Symmetrical bi-phasic square waveform when measured into a resistive load.

Environmental Specifications

Operating: Temperature Range: 0 to 35°C
Humidity: 20 to 65 % RH
Storage: Temperature Range: 0 to 55°C
Humidity: 10 to 90 % RH

Description of your unit's symbols:
There are a number of technical markings on your unit. These can be explained as follows:

Inside the battery compartment '+' indicates positive polarity and '-' indicates negative polarity. The batteries' position is also shown by 3 embossed outlines in the base of the compartment.
The unit requires three 1.5 volt (LR03 / AAA) DC batteries. DC is indicated by the symbol: 

The unit and belt are manufactured for Bio-Medical Research Ltd., Parkmore Business Park West, Galway, Ireland.
The Output Frequency indicates the number of pulses per second transmitted by the unit. This is measured in hertz ('Hz').
Output (RMSA) mean the max. output root mean square current for each channel.
Output (RMSV) mean the max. output root mean square voltage for each channel.

⚠ This symbol means 'Attention, consult the accompanying documents'.

 This symbol means type BF equipment.

 This symbol on your SLENDERTONE unit is to indicate conformity to the requirements of the Medical Device Directive (93/42/EEC). 0366 is the number of the notified body (VDE).

Inside the box, on the back of the unit is the serial number specific to this unit. The letter preceding the number indicates the year of manufacture, where 'A' denotes 1995, 'G' denotes 2001, 'H' denotes 2002 etc.

The belt's batch number is represented on the belt packaging, by the number corresponding with the  symbol.

Product Type: 551

Parameter	500Ω	1KΩ	1K5Ω
Output RMSV	5.8V	9.0V	11.1V
Output RMSA	11.7mA	9.0mA	7.4mA
Output Frequency	50-60Hz.	50-60Hz.	50-60Hz.
DC Component: Approx.	0 C	0 C	0 C
Positive Pulse Width:	250-300 μs	250-300 μs	250-300 μs
Negative Pulse Width:	250-300 μs	250-300 μs	250-300 μs
Interphase Interval:	100 μs	100 μs	100 μs

NOTE: In accordance with the laws within the member states, safety testing on the device should be carried out every two years.

NOTE: In compliance with German law, product safety testing must be conducted every two years by an authorised test service.

BELT MATERIALS

Outer material: 100% Nylon
Binding: 82% Nylon, 18% Elastane
Hook & Loop: 100% Nylon
Foam: 100% Polychloroprene

 At the end of the product lifecycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Some product materials can be re-used if you bring them to a recycling point. By re-using some parts or raw materials from used products you make an important contribution to the protection of the environment. Please contact your local authorities if you need more information about collection points in your area.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health.

English

DIARY

Fill in your diary after every session to reach your weekly target. Replacement SLENDERTONE ACTIVE diaries are available on www.slendertone.com.

Sample	4/10/06 4	/ / /	6/10/06 6	/ / /	8/10/06 6	9/10/06 4	/ / /	=	20 ★	
Day	1	2	3	4	5	6	7			
Week	1	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	2	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	3	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	4	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	5	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	6	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	7	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	8	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	9	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★
	10	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	/ / /	=	★

PRODUCT GUARANTEE

Should your unit develop a fault within two years of purchase, SLENDERTONE will undertake to replace or repair the unit or any parts found to be defective with no charge for labour or materials *, provided the unit:

- Has been used for its intended purpose and in the manner described in this instruction manual.
- Has not been connected to an unsuitable power source.
- Has not been subjected to misuse or neglect.
- Has not been modified or repaired by anyone other than an approved SLENDERTONE agent.

This guarantee complements existing national guarantee obligations and does not affect your statutory rights as a consumer.

* This excludes consumables (e.g. pads, belt etc.) when subject to normal wear and tear.

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