SLENDERTONE®

Optimum

United Kingdom: 0845 070 77 77
Republic of Ireland: 1890 92 33 88
France: 0810 347 450
Deutschland: 0800 000 1094
España: 900994467
Italy: 0800 928 388
International: +353 94 902 9936

info@slendertone.com

Designed by & Manufactured for:

Bio-Medical Research Ltd., Parkmore Business Park West, Galway, Ireland.

Copyright: © 2012, Bio-Medical Research Ltd. All Rights Reserved

Part No: 2400-5802 Rev.: 2 Date of Issue: 08/12

www.slendertone.com

SLENDERTONE®

Optimum

INSTRUCTION MANUAL

MODE D'EMPLOI

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www.slendertone.com

Contents

1.	Welcome / Product Warranty.....	3
2.	SLENDERTONE OPTIMUM kit contents.....	4
3.	How Slendertone Optimum works.....	5
4.	Contraindications and side effects.....	6
5.	Warnings.....	7
6.	Unit Controls.....	10
7.	Setup.....	11
8.	Battery charging.....	16
9.	Slendertone Optimum maintenance.....	16
10.	Frequently Asked Questions.....	17
11.	Electrostimulation and stimulation intensities.....	18
12.	Sport programmes.....	19
13.	Toning programmes.....	25
14.	Health programmes.....	29
15.	Technical specifications.....	34

Welcome

Thank you for purchasing Slendertone OPTIMUM. Please read this manual fully before using the product.

Slendertone OPTIMUM has been developed by Biomedical Research Ltd, a medical company with over 50 years of experience in Sport, Toning and Pain Relief environments.

This product uses Electronic Muscle Stimulation (EMS) technology to strengthen and rehabilitate muscles. EMS is used in medical and sport professions to improve muscle performance.

Slendertone OPTIMUM has 3 main programme groups, with 23 programmes in total, and has been designed to cater for various uses and treatments.

1. Sport... programmes build strength, speed and endurance for your chosen sport.
2. Toning... programmes shape and define your muscles creating a beautiful body.
3. Health... TENS programmes, pain relief technology which helps to reduce pain and relief from muscle soreness and daily stress.

Product Warranty

Should your unit develop a fault within two years of purchase, SLENDERTONE will undertake to replace or repair the unit or any parts found to be defective with no charge for labour or materials*, provided the unit:

- Has been used for its intended purpose and in the manner described in this user manual and accompanying Pad Placement Chart
- Has not been connected to an unsuitable power source
- Has not been subjected to misuse or neglect
- Has not been modified or repaired by anyone other than an approved SLENDERTONE agent

This guarantee complements existing national guarantee obligations and does not effect your statutory rights as a consumer.

*this excludes consumables (e.g. pads, etc.) which are subject to normal wear and tear.

Slendertone Optimum Kit Contents

Slendertone Optimum Unit

The Unit generates signals which are sent, via the pads, to your muscles to make them contract. The unit is rechargeable and takes approximately 3 hours to fully charge.

4 Dual leads

The leads are for the transmission of signals from the unit to the pads.

2 Packs of 4 small adhesive pads;

1 Pack of 4 large adhesive pads;

These self-adhesive pads are for connecting to the dual leads. The Pad Placement Chart illustrates the correct placement of the Pads for a selection of target muscles. Insert the lead plugs into the connectors at the bottom of the unit matching the colour coding on the connectors to the colours of the buttons on the unit.

Battery charger

Plug the charger into a socket and then connect your controller to the charger to charge the battery. Please ensure you only use the charger supplied by Slendertone.

Carry bag

The carry bag protects your Slendertone Optimum product between sessions.

Instruction Literature (Instruction Manual & Pad Placement Chart)


This is a detailed guide to using your Slendertone OPTIMUM product. The Pad Placement Chart illustrates the size of pads to be used and the correct placement for a selection of target muscles.

Contact

If you have any difficulty setting up your unit or have any other questions or concerns about your Slendertone Optimum unit call the SLENDERTONE Careline (see contact details on back page).

How Slendertone Optimum works

Slendertone Optimum uses **electrical muscle stimulation (EMS)** and **transcutaneous electrical nerve stimulation (TENS)** technology. This technology has been used in hospitals and by physiotherapists for over 50 years. Here is a summary of how EMS works.

Stimulation

A gentle signal is sent through the pads to your muscles.

Contraction

As the signal strength increases, the muscle contracts exactly as in normal exercise.

Relaxation

The signal then stops, allowing your muscles to relax before the next stimulation begins.

Transcutaneous Electrical Nerve Stimulation (TENS)

TENS is a pain therapy based on the stimulation of the nerve fibres, via the application of electrical stimuli to the peripheral nervous system. This works in two ways: The “pain gate” theory, which blocks the pain signals to the brain and through the increased release of endorphins, which inhibits the transmission of pain.

This technology can be used for various purposes, from pain reduction to winding-down exercises, from muscle strengthening to massage, from isometric exercises to treating muscle soreness, to the treatment of imperfections using the beauty programmes.

Slendertone Optimum is a medical unit specifically intended for use by adults with adequate knowledge of the unit acquired from reading the manual for domestic use.

How to use Slendertone Optimum

Always stimulate muscles isometrically. When stimulating the muscles of the arms or legs, bear in mind that the muscle contraction may cause involuntary limb movement, which could cause injury to you or others.

There are just three steps to setting-up a session.

- 1) Choose the programme you wish to use
- 2) Connect the pads to the unit and place them correctly on the body following your Pad Placement Chart included in your pack.
- 3) Start the programme and select the exercise intensity you feel comfortable with. Increase the intensity level as you get used to the product.

Slendertone Optimum Programmes

Sport (7 programmes) - This product will help you prepare and recover for your chosen sport, including: running, cycling, skiing, swimming and walking, by strengthening and building up the muscles.

Toning (6 programmes) - These programmes tone and shape your muscles. Apply to the desired location and select the programme that meets your current need.

Health (10 programmes) - Aches and pains can be minimised and avoided through use of these programmes. These aches and pains can be caused by sitting at a desk, standing for long periods of time, long distance travelling in a car. There are numerous programmes in this group, so you will find one to suit you and alleviate your ache and/or pain.

Slendertone recommends you consult your doctor before using Health programmes.

Contra-Indications, Warnings & Precautions

Contra-Indications, Warnings & Precautions

Your Slendertone Optimum product is suitable for use by all healthy adults; however, as with other forms of exercise, some care is needed when using it. Always follow the guidelines below and read the Instruction manual before use. The unit can deliver a strong contraction to your muscles, so although it may be used by all healthy adults, it should only be used on healthy, uninjured muscle. Some of the points below are gender specific.

Key Guidelines - very important

- This product is intended for use on healthy adult muscle only. Keep out of the reach of children.
- Never apply the Pads to any part of the head, face, neck, throat, chest area, or genital region because of the risk of stimulating inappropriate muscles and organs.
- The Slendertone Optimum product is not intended for adjunctive therapy in the treatment of medical diseases or conditions of any kind.

Warnings

- None of the programmes are designed for injured or disease-afflicted muscles. Its use on such muscles is contraindicated.
- The work imposed on the muscles by the programmes is not intended for rehabilitation or physiotherapy.
- Always ensure that two Gel Pads FROM THE SAME LEAD are used to exercise one muscle.
- Hold your unit in your hand when using it. This allows you to control your unit at all times.
- Never use the unit in ways other than those directed in this instruction manual.

CONTRA-INDICATIONS:

Do not use if:

- You have an electronic implant (e.g. cardiac pacemaker or defibrillator - as your Slendertone Optimum product may interfere with the proper functioning of the implanted stimulator) or if you suffer from any other heart problem.
- You suffer from cancer or have cancerous lesions. The effects of electronic stimulation on cancerous tissue are unknown.
- You suffer from epilepsy.
- You have recently undergone a surgical procedure.
- Following acute trauma or fracture.
- You have any other form of electronic implant. The safety has not been established for the use of electronic muscle stimulators over electronic implants.
- You are pregnant. The safety of electronic muscle stimulation (EMS) over the pregnant uterus has not been established.
- You are under medical supervision for cognitive dysfunction as you may not be able to comply with safety instructions.
- In cases of critical ischemia of the limbs.
- The unit is in close proximity (e.g., 1 yard) to shortwave or microwave therapy equipment or you are connected to high-frequency surgical equipment (e.g. diathermy), because of the risk of unit malfunction due to electronic interference.
- Wearing the unit necessitates placement over areas at which drugs/medicines are administered (short-term or long-term) by injection (e.g., hormone treatments).

PRECAUTIONS:

- Get your doctor's permission before using your Slendertone Optimum product if:
 - You wish to use your Slendertone Optimum product and have a bad back. You should also ensure the intensity is kept low.
 - You have any serious illness or injury not mentioned in this instruction manual.
 - You take insulin for diabetes.
 - You suffer from muscle or joint problems.
 - You experience any irritation, skin reaction, hypersensitivity, or other adverse reaction. You should, However, note that some reddening of the skin under the Pads is quite normal during and for a short time after an exercise session.
- After strenuous exercise or exertion, always use a lower intensity to avoid muscle fatigue.
- Contact Slendertone or an authorized dealer if your unit is not working correctly. Do not use in the meantime.
- An effective session should not cause undue discomfort.
- For first time users, muscle stimulation can be an unusual sensation. We recommend that you begin in a seated position with low stimulation intensity settings to familiarize yourself with the sensation before progressing to higher intensity settings.
- Keep your unit out of the reach of children.
- The leads and Pads must not be connected to other objects.
- Do not use your unit at the same time as any other unit which transfers an electrical current into the body (e.g., another stimulator).
- Do not touch the Pads or lead pins while the unit is switched on.
- Do not use while driving, operating machinery or cycling.
- Cease using your unit if you are feeling light headed or faint. Consult your doctor if this happens.
- Do not over exert yourself when using muscle stimulation. Any workout should be at a level comfortable for you.
- Do not place Pads over jewellery or body piercings.
- Slendertone will not accept responsibility if the guidelines and instructions supplied with this unit are not followed.
- Use of the unit is prohibited close to flammable substances or in environments with high concentrations of oxygen.
- Warning! This unit can deliver current densities in excess of 2mA/cm² when used at a high intensity with smaller pads than those supplied with this product.

Warnings

WARNINGS:

- The long-term effects of chronic electronic stimulation are unknown.
- This unit should not be used while sleeping, because you cannot react to safety information.
- You must position the Pads and operate the unit **ONLY** as indicated in this instruction manual and Pad Placement Chart.
- Never place the Pads on the front or sides of the neck, across or through the chest (i.e., one Pad on the front of the chest and one on the back), in the genital region or on the head, because of the risk of stimulating inappropriate muscles and organs.
- Avoid any recent scars, broken or inflamed skin, areas of infection or susceptibility to acne, thrombosis or other vascular problems (e.g., varicose veins), or any parts of the body where feeling is limited.
- Avoid areas of injury or restricted movement (e.g., fractures or sprains).
- Avoid placing the Pads over metal implants.
- Your Slendertone Optimum product is not intended for medical use, for the treatment of any medical condition, or for any permanent physical changes.
- Wait before using your Slendertone Optimum product until:
 - At least 6 weeks after the birth of your baby. You must consult your doctor before use.
 - One month after an IUD contraceptive unit (e.g., coil) has been fitted. You must consult your doctor before use.
 - At least 3 months after having a Caesarean section. You must consult your doctor before use.
- The heavy days of your period have finished because vigorous abdominal exercise is not recommended at this time.
- Do not use in the bath or shower, or in an environment of elevated humidity (e.g. sauna, hydrotherapy etc.).
- Caution should be used in the following:
 - When there is a tendency to bleed internally following an injury or fracture.
 - Following recent surgical procedures when electronic stimulation may disrupt the healing process
 - Over areas of skin which lack normal sensation

Adverse reactions:

- Skin irritations and burns have been reported with the use of powered muscle stimulators.
- On very rare occasions, first-time users of EMS have reported feeling light-headed or faint. We recommend that you use the product while seated until you become accustomed to the sensation.
- In rare cases, stimulation carried out in the evening can cause some people to experience difficulty in falling asleep. If this occurs, suspend the treatment and consult a doctor.


PAD PRECAUTIONS:

- To reposition Pads during a session, always stop the program currently running, reposition the Pads as directed on the Pad Placement Chart and then restart the program again.
- Only use Pads supplied by Slendertone with your Slendertone Optimum product. Any others may not be compatible with your unit and could degrade the minimum safety levels.
- The Pads are for single person use only.
- Do not plunge the Pads into water.
- Do not apply solvents of any kind to the Pads.
- Always ensure the unit is OFF before removing the Pads.
- Apply the whole surface of the Pads firmly to the skin. Do not use Pads which do not adhere properly to the skin. You can purchase new Pads at www.slendertone.com.
- If your skin is red under the Pad after a session, do not start another session on the same area until the redness has completely disappeared.


Note: If you are in any doubt about using your Slendertone Optimum product for any reason, please consult your doctor before use.

Where NEVER to place the pads


Do NOT cut the pads


On the head or face


On the neck


On both sides of the chest simultaneously


On the front and back of the thorax simultaneously.


Over skin abrasions, broken, or inflamed skin


Pads are for single person use only. Do not share them.

Wash skin thoroughly before use

Ensure skin is dry before applying the pads.

Do not use wire to remove the pads.


Unit Controls

ON/OFF button

Select programme and muscles to stimulate

Select Sport, Toning, Health programmes, male and female


Start/Stop programme, confirm selection

Increase intensity


Decrease intensity


Setup

Charger Setup

1. Select the correct plug for your region

2. Push into place until plug clicks

Selecting the language and contrast level on the display

Press the ON/OFF button  and hold it down until the language list appears on the display. Release the buttons.

Select the chosen language by using the  buttons.

Adjust the level of contrast on the display by using the  buttons.

Press  to confirm your selections.

Setting up your Slendertone Optimum


1. Identify the programme group you wish to use, Sport, Toning or Health. Then go to that section in this instruction manual to understand each of the programmes.
2. Identify the programme most suitable for your need.
3. Identify the muscle group/body area you wish to treat. (Ref: Pad Placement Chart for location of pads.)

Setup


4. Clean the skin area to be treated.

5. Connect the pads to the leads.
ALWAYS ensure that two pads FROM THE SAME COLOUR CODED LEAD are used to exercise one muscle. It is not important which pad is attached to which connector (red/black).

6. Remove the plastic covers from the pads and place on the skin as indicated in the Pad Placement Chart. Keep these covers as you will need to replace them once the session is complete.

7. Connect the leads to the controller. Take the note that the colour on the lead needs to correspond to the colour on the intensity buttons above it.


Self Adhesive Pads

The position of the Pads is shown in the Pad Placement Chart. The Pads supplied are top quality, pre-gelled and ready for use, they are particularly suitable for electrostimulation treatments. They are very flexible making it easy to treat all areas. Remove the protective plastic cover from the Pad before use, position it on the skin as shown in the Pad Placement Chart and replace the plastic cover after use. Pad duration is determined by the skin's pH value; we recommend no more than 20 applications using the same Pad. Repeated use of the same Pads can compromise stimulation safety: do not use Pads when they no longer stick to the skin; in fact it can cause skin redness that can last for many hours after stimulation. Consult a doctor if this occurs. The Pads are Latex Free.

Setup

Operating instructions (N.B. read the whole manual before using the equipment)

1. Turn Slendertone Optimum on using the  button
2. Select the programme group (Sport, Toning, Health) using the   button
3. Select the programme using the   buttons

4. Press  to select the programme
5. Select male or female with   buttons


6. Press  to select the male or female.

Setup

7. The unit displays the following message: "connect cables and electrodes and increase intensity of the utilized channels" Connect the leads to the controller. Take note that the colour on the lead needs to correspond to the colour on the intensity buttons above it.

8. Press the  button for each used channel to increase current intensity until your personal tolerance level is reached. The display shows connected leads, pads and intensity of the channels being used.


The picture above is only used as an example and may not represent exactly the S06 Endurance programme exactly.

A bar indicating the start and end of a contraction will appear on the display whilst the electronic stimulator is working. The periods of muscular contraction (bold bar) and the periods of recovery (striped bar) are graphically visualised under the bar with a cursor. This enables the user to know exactly when a contraction starts.

Setup

- At the end of the first phase, the Slendertone Optimum cancels the selected intensity, and an intermittent signal will sound. To continue with the program, press the  button for the relevant channels again until you reach the desired level of stimulation. An audio signal will sound three times when the programme has finished.
- At the end of the programme switch off the SLENDERTONE OPTIMUM and remove the leads.
- Replace the plastic covers on to the pads and place them in their original packaging.

To simultaneously increase the intensity on all channels:

 Start programme			
Increase each channel to at least intensity 1:			
 CH1 Blue	 CH2 Yellow	 CH3 Red	 CH4 White
To simultaneously increase the intensity of all channels use the following buttons:			
 Increase		 Decrease	
Skip phase: press the  button during the working programme to skip to the next phase.			
Pause/stop programme: during the working programme, press the  button once to momentarily pause the programme.			
Press the  button again to restart the programme.			
 Switch off the machine to stop the programme.			

Warning: make sure that the controller is switched off before disconnecting the pads at the end of the treatment.

Charging The Battery

The display shows when the battery is low of power. Insert the jack of the battery charger into its socket and plug the battery charger into the main power socket. Recharge for 3 hours for a complete charge.

Do not immerse the battery in water or other liquids and do not expose it to heat sources. Do not dispose of dead or defective batteries with domestic waste; dispose of in an authorised waste collection bin. Before storing Slendertone Optimum for long periods of time, we recommend that it is fully charged and stored in a dry place, checking that it is switched off.

Do not use the battery charger if:

- The plug is damaged or parts of it are broken.
- It has been exposed to water or any other type of liquid.
- The components have been damaged by a fall.

Do not immerse the battery charger in water. Use a dry cloth to clean the battery charger. Do not open the battery charger: it does not contain repairable parts.

Warning: disconnect pads during recharging.

Slendertone Optimum Maintenance

Slendertone Optimum must be switched off and the leads removed from the connectors at the end of every session. Slendertone Optimum should be kept in the carry bag supplied, together with its accessories, and carefully stored on a secure surface and protected from situations listed in the Warnings paragraph.

Cleaning: only clean the unit using a dry soft cloth

The manufacturer states that the Slendertone Optimum unit cannot be repaired by personnel external to its own company. Any work of this nature carried out by personnel not authorised by the manufacturer will be classed as tampering with the unit, freeing the manufacturer from responsibility for the warranty and hazards that the operator or user may be exposed to.

N.B. Disconnect the leads before placing the Slendertone Optimum product in the bag. If not, the leads will bend excessively close to the connectors. This may damage the leads.

Frequently asked questions

Any type of work on Slendertone Optimum must be carried out exclusively by the manufacturer or by an authorised dealer. In any event, any presumed malfunction of Slendertone Optimum must be verified before sending the unit to the manufacturer. Here below are some typical situations:

Slendertone Optimum cannot be turned on:

- Recharge the batteries and try to switch it on again; if it still does not switch on, check that the power supply socket is working.
- Check that the ON/OFF button has been pressed correctly (held down for at least one second).

Slendertone Optimum does not transmit electric impulses:


- Check that the lead jacks have been connected fully to the pads and that the plastic cover has been removed from the pad.
- Check that the leads have been connected correctly to the controller.
- Check that the pads are fully in contact with your skin.

Slendertone Optimum switches off during operation:

- Recharge the batteries for 3 hours and start the treatment again. If the problem persists contact the manufacturer.

Slendertone Optimum is switched on but nothing appears on the display:

- Open the battery compartment, lift up the battery, disconnect the connection and remove it. Inside the battery compartment you'll find a little hole with RESET writing above: insert a ball-point pen tip or similar sized object and gently press.

Disposal

To protect the environment, the unit and its accessories must be disposed of in accordance with the law in areas equipped for that purpose or as special waste.

Product subject to WEEE regulations concerning separate waste collection: when disposing of this product, please use the designated areas for disposing electronic waste or contact the manufacturer.

Electrostimulation and stimulation intensity

Electrostimulation consists of the transmission of electric microimpulses to the human body. Fields of application of electrostimulation are: pain therapy, muscle rehabilitation, athletic preparation and beauty treatments.

The stimulation intensity is displayed on the Slendertone Optimum screen for each channel in a scale progressing from 0 to 120.

Once the initial intensity has been set, Slendertone Optimum will proceed until it reaches the end of the phase; when more than one phase is present, the intensity has to be set again at the beginning of the new phase.

The types of impulses can be classified as follows:

1. **TENS impulse:** for health programmes, the intensity should be adjusted to a level between the thresholds of perception and pain. The maximum limit is reached when the muscles surrounding the area treated begin to contract. We recommend you remain below this limit.
2. **Circulation impulse:** the “pump” effect is produced by sequential contractions. The intensity must be enough to produce these contractions: the greater the contraction is, the greater the pump effect will be. But beware: there is no benefit to be gained from enduring intensity high enough to cause pain. The first electrostimulation sessions should be carried out at a low intensity to allow the body to get used to the new sensations. That way intensity can be increased gradually and will not be traumatic.
3. **Warming-up impulse:** stimulation intensity must be increased gradually to improve the state of the muscle concerned progressively. A bit like with a car: the engine needs warming-up before it is taken to full speed.
4. **Recovery impulse:** the muscle treated must visibly contract during a training impulse. The fact that the muscle tends to stiffen and increase in volume will be visible. Intensity should be increased gradually (in the first contraction) to enable you to identify the right level of stimulation comfort. Intensity can be increased up to the personal tolerance threshold during the second training contraction; this operation is then repeated during each contraction until the workload reaches the level of intensity recommended in the description of the single programmes. We recommend that you record the level of intensity reached in order to improve the level of stimulation and consequently your performance.

Electrostimulation and stimulation intensity

- 5. Massage, winding down, active recovery impulse:** intensity must be adjusted gradually to massage the muscle treated. The level of intensity should be sufficient to obtain a comfortable massage. There is no need to endure high intensity levels in this case as it is meant to be a massage, meaning that intensity can be increased gradually without excesses.
- 6. Capillarization impulse:** increase the intensity gradually to produce constant, visible stimulation of the area treated; a medium stimulation threshold is recommended, always below the pain threshold.

Sports Programmes

IMPORTANT. Stimulation intensity: the sports programmes (Resistance, Strength) are divided into warm up phases, work phases, and recovery phases. During the warm up phase (phase 1), the stimulation intensity should be adjusted to produce a slight tickling sensation in the treated muscle, this should be enough to warm up the muscle without tiring it (18-30 mA). During the work phases (the central phases of the programmes), contraction impulses are alternated with active recovery impulses. The user needs to select the two intensities separately: for the contraction, the user should see the following on the display with a bold bar; 20-30 mA for relatively unfit people, 30-50 mA for fit people, and more than 50 mA for very fit people. We recommend adjusting the intensity of the recovery phase between one contraction and another by approximately 15-20% less than the intensity of the contraction. During the recovery phase (the last phase), the intensity should be adjusted to produce a good massaging effect of the stimulated area without producing any pain (18-30 mA).

S01 • Warm Up

This programme is suitable for use before training sessions or competitions, very useful for sports involving maximum effort right from the start. Programme duration: 15 minutes/3 phases. Stimulation intensity: medium-low, increased during the last 5 minutes.

Pad Position: **see photos 01 to 20 in the Pad Placement Chart in relation to the area that you wish to stimulate.**

Sports Programmes

S02 • Resistance

This programme is designed to increase resistance to physical stress, or rather withstand intense exertion for a longer amount of time in muscle regions subjected to stimulation. Indicated for sporting disciplines involving long and intense periods of exertion.

Stimulation intensity during the contraction: start with a low level of intensity, increasing it gradually. For practised athletes the intensity used should be enough to produce visible muscle contractions.

In the event of muscle ache after stimulation, use the S07 programme (Recovery).

S02 • Resistance		Total duration: 30 minutes	
		Biceps (photo 02/14), Triceps (photo 03/15), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).	
		Abdominals (photo 01/18), Trapezius (photo 07), Large dorsal (photo 08), Gluteus (photo 17)	
		Quadriceps/Thighs (photo 10/16), Biceps femoris (photo11), Calves (photo12), Anterior Tibial (photo13).	
Phase 1	Phase 2	Phase 3	Phase 4
5 minute warm up	10 minute. 20 second recovery 10 second contraction	10 minute. 20 second recovery 10 second contraction	5 minute recovery

Sports Programmes

S03 • Strength

The Basic strength programme is used in sport to develop basic strength, which for definition is the maximum tension that a muscle can exert against constant resistance. The contractions alternate with periods of active recovery during the work phase, allowing the muscle to be trained without subjecting it to stress and improving oxygenation of the same muscle.

The following basic procedure will enable you to obtain the first results: two sessions per week (for each muscle region) for the first three weeks at medium/low intensity (20-30mA), three sessions per week for the next three weeks at high intensity (30-50mA).

Intensity must be increased gradually treatment by treatment, without overstraining the muscles. Suspend training for a few days in the event of fatigue and proceed with the "S 07" programme.

S03 • Strength		Total duration: 30 minutes	
	Biceps (photo 02/14), Triceps (photo 03/15), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/18), Trapezius (photo 07), Large dorsal (photo 08), Gluteus (photo 17)		
	Quadriceps/Thighs (photo 10/16), Biceps femoris (photo11), Calves (photo12), Anterior Tibial (photo13).		
Phase 1	Phase 2	Phase 3	Phase 4
5 minute warm up	10 minute. 13 second recovery 7 second contraction	10 minute. 14 second recovery 6 second contraction	5 minute recovery

Sports Programmes

S04 • Speed

This programme is designed to help athletes increase speed. The exercise assumes a fast pace and the contraction is short, as is the recovery. It is usually best to complete a three-week basic strength cycle of increasing intensity before using this programme. Then continue with three weeks of fast strength three times a week at high intensity (30-50mA) during the contraction.

S04 • Speed		Total duration: 26 minutes	
	Biceps (photo 02/14), Triceps (photo 03/15), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/18), Trapezius (photo 07), Large dorsal (photo 08), Gluteus (photo 17)		
	Quadriceps/Thighs (photo 10/16), Biceps femoris (photo 11), Calves (photo 12), Anterior Tibial (photo 13).		
Phase 1	Phase 2	Phase 3	Phase 4
3 minute warm up	10 minute. 14 second recovery 6 second contraction	8 minute. 10 second recovery 5 second contraction	5 minute recovery

Sports Programmes

S05 • Explosive Strength

Explosive strength programmes increase the explosive power and speed of the muscle mass, with extremely short, strengthening contractions and very long active recovery times to allow the muscle to regain strength. Complete a three-week strength cycle (S03) before using this programme. Then continue with three weeks of explosive strength twice a week. During the contraction, the intensity must be the highest that can be endured in order to obtain maximum muscle exertion whilst involving the greatest number of fibres (over 35mA).

S05 • Explosive Strength		Total duration: 28 minutes	
	<p>Biceps (photo 02/14), Triceps (photo 03/15), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).</p>		
	<p>Abdominals (photo 01/18), Trapezius (photo 07), Large dorsal (photo 08), Gluteus (photo 17)</p>		
	<p>Quadriceps/Thighs (photo 10/16), Biceps femoris (photo11), Calves (photo12), Anterior Tibial (photo13).</p>		
Phase 1	Phase 2	Phase 3	Phase 4
3 minute warm up	10 minute. 24 second recovery 6 second contraction	10 minute. 24 second recovery 6 second contraction	5 minute recovery

Sports Programmes

S06 • Endurance

The Endurance programme is used in sports to increase muscle resistance, acting mainly on slow-twitch fibres. Programme indicated for endurance sports: marathon runners, cross-country skiers, ironman, etc. Programme duration: 50 minutes. Stimulation intensity during the contraction: if not particularly fit, start with a low intensity then increase gradually. For practised athletes the intensity used should be enough to produce visible muscle contractions.

S06 • Endurance		Total duration: 50 minutes	
	Biceps (photo 02/14), Triceps (photo 03/15), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/18), Trapezius (photo 07), Large dorsal (photo 08), Gluteus (photo 17)		
	Quadriceps/Thighs (photo 10/16), Biceps femoris (photo 11), Calves (photo 12), Anterior Tibial (photo 13).		
Phase 1	Phase 2	Phase 3	Phase 4
5 minute warm up	20 minute. 20 second recovery 10 second contraction	20 minute. 22 second recovery 8 second contraction	5 minute recovery

S07 • Recovery

Can be used for all sports, after competitions or the most demanding training sessions, in particular after long and intense exertion. To be used immediately after exertion. Helps drainage and winding down, improving muscle oxygenation and helping to discharge synthetic substances produced during exertion.

S07 • Recovery		Total duration: 25 minutes	
Muscles: all the muscles you want to stimulate in a sequence			
Phase 1	Phase 2	Phase 3	Phase 4
3 minute warm up	10 minute. 6 second sequential contractions on the 4 channels	10 minute. 6 second sequential contractions on the 4 channels	2 minute recovery

Toning Programmes

E01 • Refirming

Indicated for firming up muscles; working mainly on slow-twitch fibres. Suitable for those who have never done any physical activity or have been inactive for a long period of time. Method of use:

1. Identify the muscle to be treated. To obtain good results we suggest to treat just a few muscles at a time and complete the process described below;
2. Position the pads as shown in the photo;
3. Increase the intensity until the impulse can be felt (use a low intensity for the first session to help you to understand how the machine works);
4. During the programme and over the next few days the intensity should be increased gradually so that muscle contractions are not painful;
5. During contraction generated by the unit, contract the muscle voluntarily;
6. A cycle of 15/20 applications must be completed before the first results can be seen; one application for each muscle every two days with a day of rest in between. It is possible to work on pairs of muscles, for example thighs and abdominal muscles, treating thighs one day and abdominal the next day. Working on too many muscles at the same time is not recommended.

E01 • Refirming		Total duration: 30 minutes	
	Biceps (photo 02/14), Triceps (photo 03/15), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/18), Trapezius (photo 07), Large dorsal (photo 08), Gluteus (photo 17)		
	Quadriceps/Thighs (photo 10/16), Biceps femoris (photo11), Calves (photo12), Anterior Tibial (photo13).		
Phase 1	Phase 2	Phase 3	Phase 4
5 minute warm up	10 minute. 13 second recovery 7 second contraction	10 minute. 13 second recovery 7 second contraction	5 minute recovery

Toning Programmes

E02 • Definition

Indicated for defining muscles, working mainly on explosive fibres. Suitable for people who already practise good physical activity and wish to define their muscles in greater detail. Method of use:

1. Identify the muscle to be treated. To obtain good results we suggest to treat just a few muscles at a time and complete the process described below;
2. Position the pads as shown in the photo;
3. Increase the intensity until the impulse can be felt (use a low intensity for the first session to help you to understand how the machine works);
4. During the programme and over the next few days the intensity should be increased gradually so that muscle contractions are not painful;
5. During contraction generated by the unit, contract the muscle voluntarily;
6. A cycle of 15/20 applications must be completed before the first results can be seen; one application for each muscle every two days with a day of rest in between.

It is possible to work on pairs of muscles, for example thighs and abdominal muscles, treating thighs one day and abdominal the next day. Working on too many muscles at the same time is not recommended. Stimulation intensity: 20-35mA.

E02 • Definition		Total duration: 25 minutes	
	Biceps (photo 02/14), Triceps (photo 03/15), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/18), Trapezius (photo 07), Large dorsal (photo 08), Gluteus (photo 17)		
	Quadriceps/Thighs (photo 10/16), Biceps femoris (photo11), Calves (photo12), Anterior Tibial (photo13).		
Phase 1	Phase 2	Phase 3	Phase 4
5 minute warm up	10 minute. 13 second recovery 7 second contraction	5 minute. 14 second recovery 6 second contraction	5 minute recovery

Toning Programmes

E03 • Toning

Due to a combination of capillarizing and toning impulses, this program helps mobilise fat in areas where it tends to accumulate. The Pads should be positioned to form a square around the treatment area and can be applied daily using a medium intensity.

Program duration: 30 minutes / one phase. Recommended stimulation intensity: medium.

Pads position: see photos 01 to 20 in the Pad Placement Chart.

E04 • Relax massage 1

Can be used to help relieve daily stress and after demanding training sessions To be used immediately after exertion. Helps drainage and capillarization, improving muscle oxygenation and helping to discharge synthetic substances produced during exertion. Programme duration: 15 minutes. Stimulation intensity: medium-low, increased during the last 10 minutes. Pad Position: see photos 01 to 20 in the Pad Placement Chart.

E05 • Relax Massage 2

Can be used for all sports, after competitions or the most demanding training sessions, in particular after long and intense exertion. To be used immediately after exertion. Helps drainage and capillarization, improving muscle oxygenation and helping to discharge synthetic substances produced during exertion. programme duration: 15 minutes. Stimulation intensity: medium-low, increased during the last 10 minutes. Pad Position: see photos 01 to 20 in the Pad Placement Chart.

Toning Programmes

E06 • Circulation

This programme increases microcirculation within and around the treated muscle fibres and create rhythmic contractions, facilitating the discharge of algogenic substances and promoting lymphatic activity. This can also be applied to older people to improve blood and lymphatic circulation. The programme produces sequential tonic contractions in phases 2 and 3 that reproduce the typical effect of electronic lymphatic drainage. However in this case, unlike the lipolysis programme, the contractions are sequential: first CH1 then CH2, CH3, and CH4. There are no real limits of application for this programme, which can be practiced until the desired result has been achieved. Stimulation intensity must be sufficient to produce good muscle contractions during the treatment but not enough to cause any soreness (20-30 mA).

The first results can usually be seen after 3/4 weeks of practicing 4/5 sessions a week.

E06 • Circulation		Total duration: 25 minutes	
<p>Muscles: All the muscles you want to stimulate in a sequence. Connect a channel to each muscle that you want to stimulate.</p>			
Phase 1	Phase 2	Phase 3	Phase 4
3 minute warm up	10 minute. Alternated contraction on 4 channel, 6 seconds each	10 minute. Alternated contraction on 4 channel, 6 seconds each	2 minute recovery

Health Programmes

WARNING

Pain indicates various types of medical conditions! Please consult your doctor to identify the medical condition before using Slendertone OPTIMUM.

Slendertone OPTIMUM Health Programmes deliver TENS treatment which can be used as pain therapy.

Position of pads and intensity levels

The pads should be positioned to form a square surrounding the painful area using Channel 1 and Channel 2 as shown in illustration 1 (red up/black down are not important for therapeutic purposes, follow instructions on Pad Placement Chart). The intensity should be adjusted to a level between the thresholds of perception and thresholds of pain: the maximum intensity level is the moment in which the muscles on the area treated begin to contract; over this limit, stimulation does not become more effective, just more irritating.

illustration 1

WARNING

Create a square area with the electrodes, over the painful zone. Keep 4cm minimum distance between electrodes.

Health Programmes

A01 • TENS 1

This programme also called conventional TENS, used for analgesic purposes, blocking pain at the dorsal horn of the spinal cord, in accordance with the "Gate Control Theory" by Melzack and Wall. Pain impulses leave part of the body (for example the hand) and run along the nerve tracts (through small-diameter nerve fibres) until they reach the central nervous system where the impulses are interpreted as pain. Conventional TENS activates large-diameter nerve fibres, blocking the path of small-diameter nerve fibres at the spinal cord. So action is mainly taken against the symptom: to simplify it further, the wire transmitting pain information is obstructed. Conventional TENS is a current that can be used to treat general everyday pain. The average number required to benefit from the treatment is 10/12 once a day (there are no contra indications for up to double this amount). The programme has a duration of 30 minutes in a single phase. The programme can be repeated at the end of the session for particularly persistent pain. The nature of the impulse means that the patient may experience an "addictive" effect due to which the impulse will be felt less and less: if necessary the intensity can be increased by one level to counter this effect.
Pad Position: form a square over the painful area as shown in illustration 1.

A02 • TENS 2

This type of stimulation produces two types of effects according to how the pads are positioned: positioning the pads in the dorsal region, see photo 07 in the Pad Placement Chart, promotes the endogenous production of morphine-like substances capable of raising the pain perception threshold; positioning the pads to form a square over the painful area as shown in illustration 1 produces a vascularizing effect. Vascularization increases arterial flow and consequently aids the removal of algogenic substances and helps to restore normal physiological conditions.
Treatment duration 20 minutes in a single phase, daily frequency.
Do not position the pads close to inflamed areas.
Intensity adjusted for good sollicitation of the part stimulated, the sensation must be similar to that of a massage.

Health Programmes

A03 • Cervical Pain

Specific programme for the treatment of pain in the neck area.

The intensity should be adjusted to a level between the thresholds of perception and pain: the maximum intensity level is the moment in which the muscles surrounding the area treated begin to contract; over this limit stimulation does not become more effective, just more irritating, so it is best to stop before that point. The first benefits can be seen after 10 to 12 treatments carried out on a daily basis; proceed with the treatment until the symptoms pass. Pads position: photo 21 in the Pad Placement Chart.

A03 • Cervical Pain		Total duration: 33 minutes
Pad positioning: cervical zone (photo 21)		
Phase 1	Phase 2	
TENS microimpulses 30 minutes	Decontracting 3 minutes	

A04 • Lumbar Pain

Specific programme for the treatment of pain in the lumbar area. The intensity should be adjusted to a level between the thresholds of perception and pain: the maximum intensity level is the moment in which the muscles surrounding the area treated begin to contract; over this limit stimulation does not become more effective, just more irritating, so it is best to stop before that point. The first benefits can be seen after 15 to 20 treatments carried out on a daily basis; proceed with the treatment until the symptoms pass. Programme duration: 35 minutes. Pads position: photo 23 in the Pad Placement Chart.

A04 • Lumbar Pain		Total duration: 35 minutes
Pad positioning: Lumbar area (photo 23)		
Phase 1	Phase 2	
TENS microimpulses 30 minutes	TENS burst 5 minutes	

Health Programmes

A05 • Sciatica

Specific programme for the treatment of pain in the sciatic area. The intensity should be adjusted to a level between the thresholds of perception and pain: the maximum intensity level is the moment in which the muscles surrounding the area treated begin to contract; over this limit stimulation does not become more effective, just more irritating, so it is best to stop before that point. The first benefits can be seen after 15 to 20 treatments carried out on a daily basis; proceed with the treatment until the symptoms pass. programme duration: 30 minutes. Pads position: photo 24 in the Pad Placement Chart.

A06 • TENS 3

A conventional high frequency TENS programme characterised by a narrow pulse width. The programme is suitable for acute episodes of pain. The intensity should be adjusted to a level between the thresholds of perception and pain:

Number of treatments: until pain is lessened, on a daily basis (even 2/3 times a day).

A06 • Tens 3		Total duration: 30 minutes
Pad positioning: Ankle (photo 28)		
Phase 1	Phase 2	
TENS rapid 15 minutes	Haematomas impulse 15 minutes	

A07 • Anti-Inflammatory

This programme recommended for inflammatory conditions. To be applied until the inflammatory state is lessened (10-15 applications, once a day; the daily treatments can be doubled if required). Identify the area to be treated and position the pads as shown in illustration 1(start of the Health Programme section of this manual). Adjust the intensity until a tingling feeling is produced in the area treated; avoid contracting the surrounding muscles. Programme duration: 30 minutes.

A08 • Vascularization

This programme has a vascularizing effect on the area treated. Vascularization increases arterial flow and consequently aids the removal of algogenic substances and helps to restore normal physiological conditions. Do not position the pads close to inflamed areas. Daily application is recommended, the number of applications is not defined; the programme can be used to reduce pain. Stimulation intensity should be between the thresholds of perception and slight discomfort. Programme duration: 20 minutes. Pads Position: see photo 21 and 29 in the Pad Placement Chart.

Health Programmes

A09 • Recovery

Can be used for all sports, after competitions or the most demanding training sessions, in particular after long and intense exertion. To be used immediately after exertion. Helps drainage and capillarization, improving muscle oxygenation and helping to discharge synthetic substances produced during exertion. programme duration: 35 minutes. Stimulation intensity: medium-low, increased during the last 10 minutes. Pads Position: see photos 01 to 20 in the Pad Placement Chart.

A09 • Recovery		Total duration: 35 minutes
Pad positioning: Form a square above the painful area as shown in illustration 1		
Phase 1	Phase 2	
TENS narrow pulse 15 minutes	TENS impulse length 20 minutes	

A10 • Regeneration

This programme can be used for muscle recovery and has four phases with three different impulses to vascularize and tone the superficial muscles. It also produces a local pain-relieving action. The intensity should be adjusted to a level between the thresholds of perception and pain.

A10 • Regeneration				Total duration: 40 minutes
Pad positioning: Form a square above the area to treat as shown in illustration 1				
Phase 1	Phase 2	Phase 3	Phase 4	
10 minutes Vascularization	10 minutes Tens microimpulses	10 minutes Tens sequential	10 minutes Vascularization	

Technical Specifications

Description of your unit's symbols:

There are a number of technical markings on the Slendertone Optimum unit. These can be explained as follows:

The unit requires a 1 X 4.8 volt (NiMH) DC battery pack. DC is indicated by the symbol: ---

The Output Frequency indicates the number of pulses per second transmitted by the unit. This is measured in hertz, which is indicated by "Hz".

Output (RMSA) means the maximum output root mean square current for each channel.

Output (RMSV) means the maximum output root mean square voltage for each channel.

Related outputs – Voltage/Currents: (Type 563)

Parameter	500Ω	1KΩ	1.5KΩ
Output RMSV	16.7V	33.3V	49.9V
Output RMSA	33.3mA	33.3mA	33.3mA
Output Frequency	1-110Hz	1-110Hz	1-110Hz
DC Component: Approx.	0C	0C	0C
Positive Pulse Width:	50-350 μ S	50-350 μ S	50-350 μ S
Negative Pulse Width:	50-350 μ S	50-350 μ S	50-350 μ S
Current Level:	0-120mA	0-120mA	0-120mA
Power:	0.55W	1.11W	1.66W

Risk class: IIa in accordance with the 93/42/EEC directive (and following modifications).

Insulation class: Internally powered device with type BF applied parts, in accordance with the EN 60601-1 standard.

Classification: IP30, based on the penetration of liquids and powder; equipment not protected.
Housing: ABS.

Battery charger: model Type: FW7219/NI 4-10 ,
Input 100-240 VAC / -0.3A/ 50-60Hz
Output 4.8 - 12V DC/0.8A

Self-adhesive electrodes

Limitations: equipment not suitable for use in environments with a high concentration of oxygen and/or flammable liquids and/or flammable gas; do not use with equipment for electrosurgery or short-wave or microwave therapy.

Operation: equipment suitable for continuous operation.

Output channel: 4 channels, independent and galvanically isolated.

Type of impulse: compensated biphasic square wave.

Display: reflective LCD display.

Control keypad: 10-key keypad in ABS.

Dimensions: 143X83.8X34.5 mm.

Weight: 283 gr. including batteries.

Intended Use: Muscle Stimulator

Environmental Specifications:

Operating/Storage /Transport: Temperature: from 0°C to +40°C

Humidity: 10 - 90% RH

Atmospheric pressure: from 700hPa to 1060hPa

Product expected lifetime: 5 years

This symbol on your Slendertone Optimum unit is to indicate conformity to the requirements of the Medical Device Directive (93 / 42 / EEC). 0366 is the number of the notified body (VDE).

Internally powered device with type BF applied parts

The unit is manufactured for Bio-Medical Research Ltd., Parkmore Business Park West, Galway, Ireland.

Class II equipment

Product subject to WEEE regulations concerning separate waste collection

Read the instructions for use carefully before using the device

Warning or Caution

SN stands for 'serial Number'. The Serial number is on the back of the unit. The letter preceding the number indicates the year of manufacture, where "L" denotes 2006, "M" denotes 2007 etc.

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L'apparecchio è fabbricato per Bio-Medical Research Ltd., Parkmore Business Park West, Galway, Ireland.