YU700

Tecnorita
by **BH**

Instrucciones de uso y garantía
Operating instructions and guarantee
Mode d'emploi et garantie
Instruções de utilização e garantia
Manuale d'uso e garanzia
Bedienungsanleitung und Garantie

ESPAÑOL

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www.tecnovita.eu

ENGLISH

We would like to thank you for purchasing **HairSonix** and for the trust you have placed in **Tecnovita by BH**. By means of this handy User Manual and following its simple instructions you will achieve results comparable to those obtained from the hands of a professional beautician.

We strongly recommend that you read this User Manual thoroughly. Reading and understanding it will help you use **HairSonix** effectively and safely.

Contents of box

Inside this box you will find everything you need to enjoy **HairSonix** from now onward. Carefully check through all of the items that the package should contain.

Hairsonix
Mains cable
User manual
Ultrasound gel

Basic characteristics of HairSonix

HairSonix is a Hi-Tech scalp treatment device designed for home use.

The high frequency ultrasonic wave technology helps to improve the cell membrane metabolism, contributing to better blood circulation. This technology increases the oxygenation of skin cells providing a completely safe and painless way to combat the sagging skin and it helps to increase skin metabolism.

HairSonix operating instructions

Precautions: Consult your doctor before beginning to use this appliance. This warning is especially important for people with health problems. Use by those under 18 years of age is not recommended. Keep this appliance out of children's reach.

- This appliance has been designed for home use.
- This appliance has been designed only for use on the scalp.
- **Do not use the head on your throat, neck, eyeballs, ears or ear canals.**
- Before use remove all accessories such as necklaces, earrings, contact lenses, body piercings, etc..
- Never use the appliance without first applying plenty of ultrasound gel. If you do, then the ultrasound will not be transmitted correctly and it can also damage the appliance.
- Do not stop moving the appliance. Never let it rest on the skin as this could cause burns.
- Do not use the appliance for more than 10 minutes.
- Do not use the appliance if the cable is damaged.
- Do not unplug the appliance by tugging on the cable.
- Do not move the appliance by pulling on the mains cable.
- Do not touch the appliance with wet hands
- Always clean the head after use, using a cloth slightly dampened with water. Never immerse the appliance in water.
- Do not expose the unit to high temperatures or magnetic fields.
- If you notice that the appliance is not operating correctly, do not attempt to repair, dismantle or modify it. Unauthorised repairs may cause injury, in addition to breaking the terms of the warranty. Contact your technical assistance service.

Description

- A.- Head
B – Hairsonix
E.- Adaptor

Use

1. Plug the adaptor (E) into the mains socket (220V)
2. Plug the adaptor (E) to Hairsonix.
3. Select the intensity between low (upper lights on) and high (upper and lower lights on).

Beginning a session

- Before beginning a session, remember to remove makeup and all accessories such as necklaces, earrings, contact lenses, piercings, etc. Make sure that your hair is clean.
- Always apply plenty of gel on the area and start by slowly sliding the head without pressing and keeping the same movement.

If you need to purchase ultrasound gel, you can find it in the shop where you purchased HairSonix or at our website: www.tecnovita.eu

For correct, safe use of HairSonix

For correct and safe use of HairSonix, you should keep the following application times in mind:

- **Per daily session: 10 minutes maximum.**

By session we mean the application of HairSonix on different work areas in the same day.

- **Between sessions: 48 hours minimum.**

By “between sessions” we mean the rest time between each session.

E.g.: In one week, a session schedule would be: Monday, Wednesday, Friday and Sunday.

- **Every 24 sessions (approx. 1 and a half months of treatment): rest 1 and a half months.**

After this rest period you can repeat the treatment.

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***** And, above all, respect the rest period of 48 hours between sessions. This is the time that your skin needs to rest.**

Frequently asked questions about HairSonix

¿What is HairSonix?

The high frequency ultrasonic wave technology helps to improve the cell membrane metabolism, contributing to better blood circulation. This technology increases the oxygenation of skin cells providing a completely safe and painless way to combat the sagging skin and it helps to increase skin metabolism.

What are the advantages over other methods?

5MHz ultrasonic waves are safe, painless and non-invasive and they also aid sonophoresis, allowing the nutrients in skin products to penetrate deeper.

Is it painful?

It is completely painless and safe.

Is it normal to hear a slight buzzing in your ears while using HairSonix?

Yes, some users may detect a slight buzzing sound during its application.

Can it be used on all parts of the body?

The 5Mhz HairSonix technology has been specifically designed for use on the scalp. **Under no circumstances must HairSonix be applied to the throat, neck, eyeballs, ears and ear canals.**

After how many sessions will I begin to see results?

You should notice results after the first five sessions, however results may vary depending on the skin of each individual.

How long should I wait between sessions?

It is very important to respect the 48 hour break between sessions, as this is the time that your skin needs to rest.

Are there any side effects?

If you are healthy, HairSonix has no reported side effects, nevertheless in some people it may produce a reddening of the skin that can last for a few minutes.

What is ultrasound gel for? When is it used?

The ultrasound gel is essential, ALWAYS use it before beginning a session. This will ensure that the ultrasonic waves penetrate correctly and prevent damage to the unit.

Can I use HairSonix without applying gel?

Never, because the gel acts as a conductor for the ultrasonic waves on your face. If you do not use gel then the HairSonix application loses its effectiveness and the machine can be severely damaged.

Can I use any type of gel?

We recommend using only *TECNOVITA by BH* gel, as the use of any other type may damage the machine.

Is it necessary to drink water?

No, it is not strictly necessary but it is advisable to hydrate yourself. The recommended quantity is at least 2 litres during treatment, as this helps to hydrate the skin.

When is it inadvisable to use HairSonix?

- if you suffer from hyperpigmentation, broken capillaries
- if you have sunburn
- if you suffer from severe skin problems
- active skin infections
- facial paralysis
- If you are pregnant or think you might be
- If you are breast feeding.
- undergoing fertility treatments.
- couperosis
- temporal arthritis
- malignant tumours
- diabetes
- kidney problems
- physical or mental deficiencies
- metal implants or prostheses in the area of application
- scarring in the area of application
- physical malformations.
- epilepsy
- long-term treatment with medications: anti-inflammatories, anticoagulants, antibiotics or antihistamines.
- treatment for psychological or psychiatric disorders

Can I use HairSonix if I am pregnant?

NEVER use HairSonix if you are pregnant or suspect that you might be.

Can HairSonix be used during lactation?

Do not use it if breastfeeding.

Can HairSonix be used with skin conditions?

No, It should not be used if you have any serious skin condition. Consult your doctor.

Can I use HairSonix if I have hearing implants?

No, it should not be used if you have hearing implants fitted.

Can I use HairSonix if I have implants in my mouth, gold teeth, metal implants or similar?

Do not use it if you have these types of implant, as it may cause discomfort or tenderness in the area of application.

Can I use HairSonix when I am bathing or showering?

Under no circumstances should the handset or the machine come into contact with water due to the risk of electrocution.

How do I clean HairSonix?

Never immerse HairSonix in water to clean it. Use a cloth slightly dampened with water.

Warranty

We suggest that you read through this document carefully in order to fully understand the cover that the manufacturer offers as well as your responsibilities as owner.

GUARANTEES

GENERAL INFORMATION

EXERCYCLE, S.L. would like to thank you for purchasing this product and the trust you have placed in us. Our products have been subjected to rigorous quality checks. Nevertheless, should you experience any problem with this appliance during the guarantee period then please read through these warranty conditions and then get in touch with your local authorised dealer.

This guarantee certificate has been drafted in accordance with that stated in Directive 1999/44/EC of the European Parliament and of the Council 25th May, 1999, regarding certain aspects of the sale of consumer goods and associated guarantees, incorporated into Spanish Law by virtue of Law 23/2003, 10th July; Guarantees on the sale of consumer goods.

To this effect, EXERCYCLE, S.L. VAT nº B-01029032, with business address c/ Zurrupitieta 22, Vitoria, Spain, guarantees this product under the following terms and conditions:

1. GUARANTEE PERIOD

The product covered by this guarantee certificate is protected for a period of TWO (2) YEARS as from the date of receiving the product, which will be taken as the date that appears on the sales receipt for the product or on the corresponding delivery docket should this be later.

It is essential that this be presented in the event of possible repairs. The policy regarding the warranty period may vary between different countries.

The stated warranty period is only applicable to Europe.

2. RIGHTS UNDER THE GUARANTEE

- Repair or replacement of the product which, in such instance, will be carried out within a reasonable period and without significant inconvenience to the consumer, taking into account of the nature of the goods and the purpose for which the consumer required them.
- Those cases in which the repair or replacement represents a disproportionate cost for the seller in comparison with other forms of remedy will be exempt assuming these do not prejudice the consumer.
- Repairs will be free for the consumer with regard to labour, materials and carriage.
- In the event that the machine is repaired or replaced with a new one, the guarantee period will be suspended temporarily until the new appliance has been delivered and will carry at least six months of guarantee, however under no circumstances will a new 2 year period be initiated.

The buyer may, as is apparent, require (i) a reduction in price or proceed to a (ii) recession of the contract when the repair or replacement is not carried out within a reasonable time or with significant inconvenience to the buyer. The reduction in price (i) will be proportional to the difference between the value that the goods would have had if they had been delivered in accordance with the contract and the actual value that the goods were worth at the time of the said delivery. Recession of the contract (ii) will not proceed when the non-conformity is of minor importance.

3. THE GUARANTEE DOES NOT COVER

- Defects attributable to use other than that established in the User Manual.
- Professional use or commercial leasing or hiring of the product.
- Poor installation, poor assembly or breakage of components due to improper installation, (the User Manual gives detailed instructions on how to install and use the machine).
- Improper use of the product (use without applying gel, exceeding operating times, exposure to high moisture levels, etc.) or problems arising from lack of maintenance (accumulation of dirt on any part of the product or its accessories, etc.). (The User Manual provides advice on how to install and use the product).
- Repairs carried out by persons not authorised by EXERCYCLE, S.L. or appliances that have been opened or tampered with by persons not authorised by the manufacturer.
- Deviations from the standard product, such as light vibrations, noises or breakage of components that do not affect the performance of the machine and which cannot be attributed to manufacturing faults.
- Scratches, scrapes, abrasions, knocks, falls, etc. caused by carelessness.
- Defects on the appliance caused by damage during transit.

4. CLAIMING AGAINST THE GUARANTEE

- Claims made for nonconformities will be limited to three years as from the date that the goods were delivered.
 - The consumer must report the existence of the defect within a period of 2 months as from the date that the defect was discovered.
 - This guarantee does not affect the consumer's statutory rights under any applicable law, unless it is incompatible with the application of other actions derived from compensation sought for hidden faults in the contract of sale.
 - With regard to queries that may arise from the interpretation or application of this guarantee, these will be governed by the corresponding national law.
- Any discrepancy arising from the interpretation, validity or application of this contract will be subject to the courts and tribunals of Vitoria, Spain.

Certificates

By virtue of this document, Exercycle S.L. declares that this product meets the essential requirements and other relevant provisions of Directives 2004/108/EC and 2006/95/EC

Technical specifications

Weight: 0.2 kg

Dimensions: L x An x Al: 13 x 8 x 5 cm

Ultrasound: 5 Mhz

Voltage: AC 220 V

Output: 9V – 1000mA

Please keep the receipt as proof of purchase.

Please keep the receipt as proof of purchase.

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