Manufacturer:
DJO, LLC
1430 Decision Street
Vista, CA 92081 U.S.A.

EU Authorised Representative
MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany

© 2010 DJO, LLC

One


Designed by Compex®

QUICK START GUIDE


Description of stimulator

- A** On/Off button
- B** Sockets for 4 stimulation cables
- C** Stimulation cables
Channel 1 = blue Channel 2 = green Channel 3 = yellow Channel 4 = red
- D** +/- buttons for 4 stimulation channels
- E** Socket for the charger (slide the red cover to the right to free the charger connector)
- F** Battery compartment
- G** Belt clip socket

1. Connecting the Cables


2. Selecting Language, Contrast and Volume

- B** Use the channel 1 +/- button to choose the language you wish to use.
- C** Use the channel 2 +/- button to adjust the contrast of the screen.
- D** Use the channel 3 +/- button to adjust the volume.
- E** Use the channel 4 +/- button to adjust the backlight.
On: Backlight always on.
Off: Backlight always off.
Auto: Backlight activated whenever a button is pressed.
- A** Use the On/Off button to confirm and save your choices. The settings will be applied immediately.


3. Selecting a category

- A** Press the On/Off button to turn off the unit.
- B / D** Pressing +/- on channels 1 and 3 enables a category to be selected.
- E** Press the channel 4 +/- button to confirm your choice.


4. Selecting a programme

- A** Press the On/Off button to return to the previous screen.
- B** Use the channel 1 +/- button to choose a programme.
- E** Press the channel 4 +/- button to validate your selection and start the stimulation session.


5. Starting a programme

- A** Press the On/Off button to place the unit in Pause mode.
- B C D E** The device beeps and the symbols of the four channels flash, changing from + to 000: the four channels are at 000 energy. You must increase the stimulation energy so that the stimulation can start. To do this, press the + buttons for the relevant channels until the desired setting is reached.

6. Ending a programme

- A** To turn off the stimulator, press the ON/OFF button.

For a more in-depth description of the sequences introduced in this quick start guide, refer directly to your user manual.

For explanation on how to use the programmes and example of specific applications, refer directly to your user manual (chapter 7).

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